# Basic


## international: **EASyS Advantage S**

Ⓓ	<b>Bedienungsanleitung</b> – Basic Sitzeinheit .....	6
Ⓖ	<b>User manual</b> – Basic seat unit .....	14

(1)

(2)

(3)

(4)

(5)

(6)

(7)

(8)

(9)

(10)

(11)

(12)

(13)

(14)

(15)

(16)


(17)

(18)

(19)

(20)


(21)

(22)

(23)

(24)

(25)

(26)

(27)

(28)

(29)

(30)

(31)

(32)

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**Dear user of the Basic:**

The Basic seat unit you have purchased is an innovative product made by Thomashilfen.

Please read this user manual in order to ensure that the handling of the Basic seat unit – on the chassis of your choice – is safe, practical and comfortable for you, and especially for your child.

If you still have any questions or problems, please call your local dealer or contact us directly.

We want to take this opportunity to draw your attention to our website [www.thomashilfen.com](http://www.thomashilfen.com), where you can find the latest information and important product safety information, as well as information about any product recalls.

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**Table of Contents**

	<b>Page</b>
Purpose and areas of application.....	14
Safety checks and maintenance interval .....	15
General notes regarding safety .....	15
Application risks and contraindications .....	15
Symbols and warning signs .....	16
Specifications and versions .....	16
Scope of delivery .....	17
Basic equipment of the seat unit.....	17
Setup and assembly .....	17
Possible settings .....	18
Attaching the belts .....	19
Transport .....	19
Product / serial number.....	20
Cleaning, disinfection and care .....	20
Warranty terms / CE label .....	20
Maintenance schedule – annual inspection .....	21

**Purpose and areas of application**

When combined with a matching chassis, the Basic seat unit is suitable for the transport (no more than walking speed) of children with disabilities by their parents or other suitable assistants indoors and outdoors. The seat unit, in conjunction with the A chassis Advantage, must not be used for other purposes (for example, as a transport vehicle for more children, for shopping bags, etc.).

When used with the “Crash” version of the A chassis Advantage, the Basic has been approved for transporting children in motor vehicles and vans for disabled people in accordance with ISO 7176-19.

## Safety checks and maintenance intervals

A visual inspection of all the operating controls and fastening screws should be carried out each month by a qualified specialist.

Please remember to tighten screws, nuts and clamping levers after each adjustment.

### Please observe the following notes regarding safety

- Before using the product for the first time, read the user manual carefully; if you have poor eyesight, have someone read it to you. If the user manual is lost, you can download it at any time in the download area of our website [www.thomashilfen.com](http://www.thomashilfen.com). The stored PDF document can also be displayed at a larger size.
- We recommend that you adapt the settings of the Basic to your child's current body measurements at least every 3 months. If necessary, this can also be done with the help of your therapist / orthopedic technician.
- **IMPORTANT:** Always secure the child using the 2-point pelvis belt or a different retaining belt (sold separately); never leave the child unattended in the rehab stroller.
- **CAUTION:** The gap between the chassis and the brake release rod may pose a risk of getting clothes and limbs caught in the area of the feet.
- **CAUTION:** The gap between the seat and the chassis may pose a risk of getting clothes and limbs caught in the seat area.
- **CAUTION:** Sunlight, electric heaters or other heat sources may excessively warm up the seat / rehab stroller. In this case, there is the risk of burning your back and legs upon contact! Therefore, always ensure that the seat / rehab stroller is not exposed to these heat sources during use and storage; use it only when it has cooled sufficiently.
- **CAUTION:** When stopping for a long period, you should lock the rehab stroller using the foot brake, thus preventing the stroller from unintentionally rolling away.
- Please do not overload your rehab stroller and observe the maximum load capacity (see the specifications).
- In the dark, reflective clothing makes you and your child more visible for other road users.
- Attaching heavy bags or string bags to the stroller increases the risk of tipping. Use the basket on the chassis.
- When handling the Basic seat unit (especially when folding and unfolding it), heed the recommendations in the user manual.
- Please keep your children away from the packaging; there is a danger of suffocation.
- Only one child at a time may be placed in the rehab stroller.

Previous market experience has not shown any dangerous combinations of devices with the seat unit.

## Application risks and contraindications

When used appropriately, risks for the application can be excluded. Contraindications are not known. The cushioning material used in the Basic is "flame retardant" and has been tested according to DIN EN 1021-2.

The degree of resistance to inflammation in materials and components was tested according to standards and deemed to be safe.

## Symbols and warning signs

Sticker	Meaning
	<p><b>Specification labels / loading capacity stickers</b></p> <p>The specification label of the seat unit is attached at the bottom of the seat.</p> <p>The specification label includes the name of the manufacturer, including the address; the type designation, the maximum loading capacity and the CE label.</p>
	<p><b>Bar code label</b></p> <p>The bar code label is attached under the seat of the seat unit.</p> <p>The bar code label includes the item number, designation, serial number and production date of the product.</p>

## Specifications and versions

Basic seat system	Basic size 1	Basic size 2
Seat depth	190 - 300 mm / 7.5 - 11.8"	270 - 380 mm / 10.6 - 15"
Seat width	190 - 310 mm / 7.5 - 12.2"	250 - 370 mm / 9.8 - 14.6"
Back height	415 - 680 mm / 16.3 - 26.8"	550 - 740 mm / 21.7 - 29.1"
Lower leg length	175 - 290 mm / 6.9 - 11.4"	210 - 370 mm / 8.3 - 14.6"
Hip angle	90° to 130°	90° to 130°
Knee angle	90° to 180°	90° to 180°
Foot angle	78° to 100°	78° to 100°
Turning range of seat	180°	180°
Seat tilt-in-space (on A chassis Advantage) driving direction	-10° to +35°	-10° to +35°
Seat tilt-in-space (on A chassis Advantage) face-to-face	0° to +35°	0° to +35°
Size of hip guide (w x h)	230 x 120 mm / 9.1 x 4.7"	300 x 150 mm / 11.8 x 5.9"
Size of footrest (w x d)	280 x 180 mm / 11 x 7.1"	320 x 210 mm / 12.6 x 8.3"
Overall dimensions (on A chassis Advantage) (l x w x h)	995 x 595 x 1140 mm 39.2 x 23.4 x 44.9"	1070 x 690 x 1200 mm 42.1 x 27.2 x 47.2"
Folded dimensions (on A chassis Advantage) (l x w x h)	800 x 595 x 670 mm 31.5 x 23.4 x 26.4"	800 x 690 x 630 mm 31.5 x 27.2 x 24.8"
Weight (seat)	8.6 kg / 19 lb	10.2 kg / 22.5 lb
Max. user weight (seat)	35 kg / 77.2 lb	40 kg / 88.2 lb

### The Basic seat unit is available in the following versions:

Item No.	Designation	Equipment
6226	Basic seat unit size 1	Seat unit including cushioning and footrest
6426	Basic seat unit size 2	Seat unit including cushioning and footrest

## Scope of delivery

The Basic seat unit is delivered pre-assembled. If the seat unit was ordered together with the A chassis Advantage, the seat unit is lying upside down on the chassis upon delivery. The seat unit has been upholstered at the factory with seat and back cushions. The additional pad is attached under the seat cushion. In addition, the seat unit is supplied with a user manual and an Allen wrench. The Allen wrench is required for adjusting the height of the footrest.

## Basic equipment of the seat unit

Physiological pivot points in the hip and knee joints; breathable seat and back cushions as well as side pads with safety reflectors; additional pad for more comfort; body side guide; angle-adjustable, foldable and height-adjustable foot support; holder for sun and rain cover without leg blanket; infinite adjustment of the back height, hip angle, seat depth and seat width; easily accessible setting by means of star screws.

## Setup and assembly

The settings / adjustments must always be made by trained / instructed professionals (such as medical product consultants of specialist medical dealers). Using the release lever, the seat can be removed from the A chassis Advantage without tools. It is important that the release lever remains actuated while inserting and removing the seat unit. No tools are required when using the seat unit and the chassis for the first time, except for the Allen wrench that is included in the scope of delivery.

The accessories are attached / removed using commercially available tools. Separate installation and operating instructions are included with the accessories. If these are missing, you can download the instructions from the download section of our website [www.thomashilfen.com](http://www.thomashilfen.com) at any time.

### Inserting the seat unit

Activate the release lever on the backrest and insert the seat unit into the chassis from above (Fig. 1). Make sure that the buttons on the sides of the seat unit between the folding joint and the connecting strap of the rear frame and the sliding support are properly threaded through (Fig. 2). Now guide the buttons on the sides of the seat unit into the receiving slots on the sides of the chassis (Fig. 3) and let the seat unit slide into position (Fig. 4). When you release the trigger, the seat unit engages into the closest position of the hole matrix.

**CAUTION:** Always ensure that the locking pins on both sides of the seat unit engage firmly into a hole in the hole matrix.

### Removing the seat unit

To remove the seat unit, actuate the release lever and pull the seat unit upwards out of the frame (Fig. 1). Make sure that the buttons on the sides of the seat unit between the folding joint and the connecting strap of the rear frame and the sliding support are properly threaded through (Fig. 2).

### Attaching the footrest

The footrest must be fastened to the knee joints. To do this, align the ratchet joints and fasten them to the screws (from the inside), the saddles (from the outside) and the eccentric levers (from the outside; Fig. 5).

### Seat cushion

The seat cushion consists of two elements: the main gray cushion as well as a black additional pad for more comfort. The additional pad is attached with hook and loop to the underside of the main cushion (Fig. 6). In this regard, pay attention that the notch in the additional pad is facing forward (Fig. 7). Fasten the seat cover by firmly attaching it to the rear end of the seat plate using the hook and loop flap (Fig. 8). The lower two hook and loop straps are wrapped around the tubes of the footrest.

If necessary, the additional pad can be removed.

## Backrest upholstery

The backrest upholstery is attached at the top of the back plate, as is the hook and loop strip on the back plate (Fig. 9). The lower part of the backrest upholstery is inserted to the rear between the seat and the back (Fig. 10). The lower end is provided with hook and loop and is attached to the rear end of the seat plate to close the gap between the seat and the back (Fig. 11).

## Side cushion

Place each of the shorter cushions of the two side cushion pairs over the side tubes so that the hook and loop is facing inward (Fig. 12). The hook and loop fastener is located on the outside of the back plate curve. Attach the cushion in such a manner that the upper edge of the cushion on the back ends just below the thumbscrew (Fig. 13). Then attach the rear of the cushion with hook and loop (Fig. 14). The zippers of the shorter side cushion can be positioned inwards if necessary.

Hook and loop each of the longer cushions to the inside of the back plate curve. The longer cushions are attached at a height where they are flush with the top end of the outer side cushions (Fig. 15).

## Lateral seat upholstery

The lateral seat upholstery can be opened with a zipper. It is pulled over the seat side guides in such a manner that the reflective strip on the outside runs from the rear / top to the front / bottom (Fig. 16).

## Holders for sun and rain cover without leg blanket

The Basic comes with holders for the sun and rain cover (without leg blanket) as standard (Fig. 17). Their height can be adjusted as required. To do so, loosen the screws of the holders for the sun and rain cover and move it on the tubes into the desired position. Then tighten the screws firmly again.

## Possible settings

### Seat width

To adjust the seat width, loosen the four soft-touch nuts under the seat of the Basic and move the side elements accordingly (Fig. 18). Then retighten the nuts. If the seat depth is set to the minimum depth, the maximum adjustable seat width depends on the gap before the lateral seat upholstery of the back.

### Seat depth

Loosen the two front soft-touch nuts under the seat of the Basic and pull the seat plate forwards (Fig. 18). Then retighten the nuts.

**TIP:** To easily adjust the seat depth, use one hand to grip the front edge of the seat and use the other to hold the lower end of the back (in the center) steady. The depth of the seat side guide can also be adjusted. To do so, open the cushions of the seat side guides with the zippers and remove them. Loosen the two screws, set the desired depth and retighten the screws (Fig. 19).

### Back height

Loosen the two thumbscrews on the side tubes (Fig. 20). Now the height of the backrest can be infinitely adjusted up to the engagement point (spring pins engage with a click). To reduce the back height, press the spring pins below the thumbscrews and slide the back plate downwards to the desired height (Fig. 21).

To achieve the maximum back height, also loosen both thumbscrews on the upper cross bar and pull the back plate upwards (Fig. 22). Then retighten all the thumbscrews.

**TIP:** The extended tubular bar provides additional safety for children who are predominantly transported in the horizontal position.

## Angle of recline

Loosen the clamping lever on the backrest (Fig. 23) and remove the safety pin below the lever (Fig. 24). Now the backrest can be brought into the desired position. Finally, retighten the clamping lever after the position has been set.

When transporting the rehab stroller in a vehicle, the safety pin must always engage into the hole (90° seat angle).

## Lower leg length

The grub screws on the side of the footrest can be loosened using the Allen wrench (SW 3; Fig. 25). Now the footrest can be adjusted infinitely. Then the grub screws must be retightened.

## Foot angle

Loosen the four screws on the right and left of the underside of the footrest, adjust the desired foot angle, and then retighten the screws. You will need a 10 mm box wrench or combination wrench (Fig. 26).

## Knee angle

Set the knee angle by loosening the red eccentric lever on both knee joints (Fig. 5).

**TIP:** The footrest can be folded up to simplify getting into and out of the stroller.

## Attaching the belts

Depending on the size of the child, guide the two upper belts through the openings in the cushioning and insert them through the slots in the frontmost back plate (Fig. 27). Then guide the belts into the buckles and adjust the desired belt length (Fig. 28).

The two bottom ends of the belts can be attached with buckles to the rear of the back plate or to the plastic element on the side (Figs. 29-30). Guide the belts to the left and right of the backrest upholstery and continue to the corresponding buckle, and then fasten them.

The front belts of the 4-point lap belt are also attached to the plastic element on the right and left sides of the seat side guide. Remove the cushion of the seat side guide, guide the belt through the slot in the seat side guide, insert it into the buckle, and fasten it. Attach the cushion again. (Figs. 31-32)

**TIP:** In addition, the belt length can be adjusted at the front; to do so, pull on the D-ring to shorten the belt.

## Transport

The Basic seat unit, together with the A chassis Advantage, can be folded up for transport in a vehicle trunk. Heed the instructions for folding the rehab wagon and use the lock for folding to avoid unintentional unfolding.


The seat unit and chassis can be separated so that their weights can be lifted separately. In this regard, please refer to the instructions for removing the seat unit.

To ensure the simplest possible handling, use one hand to grab the seat unit by the handle on the backrest, to which the release lever for the seat tilting is also attached (Fig. 1). With the other hand, grab the front edge of the seat.

## Product / serial number

To obtain the correct spare parts or for any technical call-backs that may be required, we always need you to provide us with the serial number and the precise item number.

The serial number and the item number of the Basic seat unit are located on the bar code label that is attached under the seat.

## Cleaning, disinfection and care

### Disinfection

Before disinfection, clean all the cushions and handles. Wipe all parts of the product with a cloth that is damp with disinfectant. Only water-based disinfectants should be used. When disinfecting, the manufacturer's application instructions must be heeded.

### Cushions

All cushions can be easily removed in a few simple steps. The 100% Trevira CS material is a high-quality fabric that is washable up to 40 °C. The cushions should be washed separately; the hook and loop straps must be covered before washing. Please also refer to the care instructions.

We would like to point out that even high-quality fabrics can fade as a result of continuous and intensive sunlight or frequent washing. The cushions are flame-retardant (EN 1021-2).

### Hook and loop fasteners

In order to maintain the functionality of the hook and loop fasteners, brush off the hook and loop strips from time to time. When washing the hook and loop straps, always close and cover them.

## Warranty terms / CE label

Thomashilfen offers a 3-year warranty on all frame parts from the date of purchase. The warranty covers any claims affecting the functioning. Excluded is damage that arises due to improper handling (e.g. overloading) as well as normal wear.

The Basic meets the requirements of European standard EN 12182 and EN 12183 and is provided with a CE label.

We hope you enjoy your Basic rehab stroller!

## Maintenance schedule – annual inspection / reutilization

### Preparation

A user manual must accompany the product. If none is present, please order one from the manufacturer. Familiarize yourself with the functions of the product. If you are not familiar with the product, read the instructions carefully before the inspection. Clean the product before the inspection. If necessary, heed the care instructions and product-specific inspection instructions in the user manual.

Item	Area ☑ = Inspection performed	(I = Inspection / R = Reutilization)							
		Function		Damage		Warping		Replace	
1	Seat	I/R		I/R		I/R			
2	Angle adjustment	I/R		I/R		I/R			
3	Height adjustment	I/R		I/R		I/R			
4	Seat cushion			I				R	
5	Missing caps and plugs							I/R	
6	Back	I/R		I/R		I/R			
7	Angle adjustment	I/R		I/R		I/R			
8	Height adjustment	I/R		I/R		I/R			
9	Back cushion			I				R	
10	Push bar	I/R		I/R		I/R			
11	Push bar cover			I				R	
12	Chassis	I/R		I/R		I/R			
13	Tires	I		I		I		R	
14	Hand brake	I/R		I/R		I/R			
	Options								
15	Headrest	I/R		I/R		I/R			
16	Headrest cushion			I				R	
17	Lateral trunk support	I/R		I/R		I/R			
18	Lateral trunk support cushions			I				R	
19	Armrests	I/R		I/R		I/R			
20	Arm supports			I				R	
21	Lateral hip support	I/R		I/R		I/R			
22	Lateral hip support cushions			I				R	
23	Abduction wedge	I/R		I/R		I/R			
24	Abduction wedge cushion			I				R	
25	Lateral tibia support	I/R		I/R		I/R			
26	Lateral tibia support cushion			I				R	
27	Tray	I/R		I/R		I/R			
28	Belts	I/R		I/R		I/R			
Maintenance was carried out		on:							

Service and repairs to the rehab stroller may be carried out only by qualified personnel. The user of this rehab aid must ensure that the required inspections / maintenance be carried out regularly and on schedule. If damage is evident, the user must actively inform the qualified personnel and have the damage repaired immediately. Only original spare parts may be used for repairs (the list of replacement / spare parts can be found in the download area on our website [www.thomashilfen.com](http://www.thomashilfen.com)).

Removal / attachment of spare parts must always be performed by qualified personnel only! If it is necessary to return the product to your dealer, please make sure that you pack the chassis so that it can be transported without damage.

Member of the International  
Support Association for the  
Rehabilitation of Children and Youths e.V.

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Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG