

# THOMASHILFEN Major Order form



**Rehab Safety Car Seat for children, Allrounder with Seatfix as integrated equipment**

**Restraint system for special needs „S“  
ECE R 44/04 - Group I - III  
from 9 kg / 19.8 lbs  
to 36 kg / 79.4 lbs**

**Order e-mail:**  
info.export@thomashilfen.de

**Status: 01.01.2024**  
valid until further notice

**Address** Customer no.:

Company: \_\_\_\_\_

Street: \_\_\_\_\_

Postal Code: \_\_\_\_\_

City: \_\_\_\_\_

Commission: \_\_\_\_\_











Contact person: \_\_\_\_\_

Phone: \_\_\_\_\_

E-mail: \_\_\_\_\_

Date: \_\_\_\_\_

BASIC MODEL	item code
-------------	-----------

<input type="checkbox"/>		Major Car Seat, grey		980000610000
<input type="checkbox"/>		Major Car Seat, red		9800006110000
<input type="checkbox"/>		Major Car Seat, blue		9800006120000
<input type="checkbox"/>		Major Car Seat, green		9800006130000
<input type="checkbox"/>		Major Car Seat, magenta		9800006140000

**Basic configuration:**

Seat and back with integrated 7-fold height-adjustable headrest; seat cushion; back cushion; head cushion; additional head cushions (left, right); crotch pad, small; 5-point belt with central adjuster; chest clip; shoulder protector; cup holder (can be used right or left); removable covers, machine washable up to 40° C; Seatfix connectors for anchoring the seat in appropriately equipped vehicles

**Measurements:** seat depth: 32 cm / 12.6" // seat width (with art. 674 hip supports): front: 31,5 cm (21,5 cm) / 12.4" (8.5") // back: 26 cm (19 cm) / 10.2" (7.5") // back height (seat): 62 – 74 cm / 24.4 - 29.1"

 please find accessories on page 2 and 3

Customer no.:

Company

Contact person:

**Status: 01.01.2024**  
valid until further notice

**Order e-mail:**  
info.export@thomashilfen.de

**ACCESSORIES – can ONLY be used WITH adapter plate** item code

**BASE – An adapter plate is required to use the + accessories.**

- Adapter plate** 980000500000
- ↓  
can be combined with:
- +  **Swivel base\*** 9800005010000
- +  **Back stabilizer** (required when using swivel base) 9700663080000
- +  **Cushion for swivel base** 9800005130000
- +  **Seatfix adapter\*\*\*** 9800005060000
- +  **Stabilizing bow, back\*\*\*** 9800005070000

**Recommendation:** For a better support always use the swivel base with Seatfix adapter or stabilization bow.


- +  **Tilting base** (5,5° / 11° / 16°) 9800005020000
- +  **Seatfix adapter\*\*\*** 9800005060000
- +  **Stabilizing bow, back\*\*\*** 9800005070000

**Recommendation:** For a better support always use the tilting base with stand or footrest.

- +  **Footrest, size 1\*\*** (depth: 12 cm / 4.7"), tiltable and adjustable in angle 9800005030000
- +  **Footrest pad** (10 cm / 3.9"), can be combined with art. 507 5005110501010
- +  **Footrest, size 2\*\*** (depth: 18 cm / 7.1"), tiltable and adjustable in angle 9800005040000
- +  **Stand\*\*** 9800005050000
- +  **Seat depth extension** (+ 13 cm / + 5.1") 9800005080000
- +  **Abduction block, adjustable** (depth: +15 cm / +5.9") 9800005090000

**NOTE:**  
 \* Swivel base cannot be used in car seat group I (9 - 18 kg / 19.8 - 39.7 lb)  
 \*\* Footrest and stand cannot be used in combination  
 \*\*\* Seatfix adapter and stabilizing bow (back) cannot be used in combination  
 \*\*\*\* Seat wedge (below) and swivel base cannot be used in combination

↪ please find further accessories on page 3



# THOMASHILFEN Major Order form (accessories)

Customer no.:

Company

Contact person:

**Status: 01.01.2024**  
valid until further notice

**Order e-mail:**  
info.export@thomashilfen.de

further ACCESSORIES – can be used WITHOUT adapter plate	item code
<input type="checkbox"/> Abduction block, padded	980000670000
<input type="checkbox"/> Seat wedge, below**** (+10°)	9800006710000
<input type="checkbox"/> Seat wedge, inside	9800006720000
<input type="checkbox"/> Table	9800006730000
<input type="checkbox"/> Lateral trunk pads (foam core), grey	● 9800006750000
<input type="checkbox"/> Lateral trunk pads (foam core), red	● 9800006751000
<input type="checkbox"/> Lateral trunk pads (foam core), blue	● 9800006752000
<input type="checkbox"/> Lateral trunk pads (foam core), green	● 9800006753000
<input type="checkbox"/> Lateral trunk pads (foam core), magenta	● 9800006754000
<input type="checkbox"/> Hip supports, grey	● 9800006740000
<input type="checkbox"/> Hip supports, red	● 9800006741000
<input type="checkbox"/> Hip supports, blue	● 9800006742000
<input type="checkbox"/> Hip supports, green	● 9800006743000
<input type="checkbox"/> Hip supports, magenta	● 9800006744000
<input type="checkbox"/> Crotch pad, large, grey	● 9800006760000
<input type="checkbox"/> Crotch pad, large, red	● 9800006761000
<input type="checkbox"/> Crotch pad, large, blue	● 9800006762000
<input type="checkbox"/> Crotch pad, large, green	● 9800006763000
<input type="checkbox"/> Crotch pad, large, magenta	● 9800006764000
<input type="checkbox"/> NapUp grey	9800806550000
<input type="checkbox"/> Anti-escape vest	9800008260000
<input type="checkbox"/> Car seat protection mat incl. back	9800008141000

**NOTE:**

- \* Swivel base cannot be used in car seat group I (9 - 18 kg / 19.8 - 39.7 lb)
- \*\* Footrest and stand cannot be used in combination
- \*\*\* Seatfix adapter and stabilizing bow (back) cannot be used in combination
- \*\*\*\* Seat wedge (below) and swivel base cannot be used in combination



Abduction block, padded




Seat wedge, below



Seat wedge, inside



Table



Lateral trunk pads (foam core)  
(coloured, matching the basic model)



Hip supports (coloured, matching the basic model)



Crotch pad, large (coloured, matching the basic model)



NapUp grey



Anti-escape vest



Car seat protection mat incl. back