

RECARO Monza Nova 2 Reha Order form



Rehab Safety Car Seat for children

Restraint system for special needs „S“

ECE R 44/04 - Group II - III

from 15 kg / 33.1 lbs tested up to 50 kg / 110.2 lbs

Order e-mail:
 info.export@thomashilfen.de

Status: 05.03.2024
 valid until further notice

Address **Customer no.:**

Company: _____

Street: _____

Postal Code: _____

City: _____

Commission: _____







Contact person: _____

Phone: _____

E-mail: _____


Date: _____

BASIC MODEL	item code
-------------	-----------

<input type="checkbox"/>		RECARO Monza Nova 2 Seatfix Reha, petrol		980008015000
<input type="checkbox"/>		RECARO Monza Nova 2 Seatfix Reha, frozen blue		980008016000
<input type="checkbox"/>		RECARO Monza Nova 2 Seatfix Reha, silent grey		980008017000

Standard features: Seat and back with integrated 11-fold height adjustable headrest; including integrated Seatfix locking arms which can be fixed to the location points of vehicles which are equipped with Isofix; stepless neck support of the headrest by an integrated air cushion; side impact protection by energy absorbing elements; removable head, lateral trunk and side supports; fixed armrests; 5-point positioning belt; crotch pad; removable abduction block; integrated speakers in the headrest for mp3-player use, integrated mp3-player pocket; RECARO Air Circulation System for optimum air circulation; ventilated foam; removable covers, can be machine-washed at 30°C (86°F)

Measurements: seat depth: 30 cm / 11.8" (43 cm / 16.9" with seat depth extension) / seat width: front 29 cm / 11.4", back 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside) / back height: 60 - 77 cm / 23.6 - 30.3"

 please find further accessories on page 2

RECARO Monza Nova 2 Reha Order form

Customer no.:

Company

Contact person:

Status: 05.03.2024
valid until further notice

Order e-mail:
info.export@thomashilfen.de

ACCESSORIES	item code
-------------	-----------

Adapter is required to use the + accessories.

- | | |
|---|----------------------|
| <input type="checkbox"/> Adapter for accessories (seat depth extension, footrest, stand, stabilising bow, seatfix) (incl. adapter plate) | 9800008212000 |
| + <input type="checkbox"/> Seatfix connection ** (for swivel base) option 1 (if vehicle is equipped with Isofix) | 9800008181000 |
| + <input type="checkbox"/> Stabilising bow, rear (for swivel base) option 2 (if the vehicle is NOT equipped with Isofix) | 9800008190000 |
| ↓ can be combined with: | |
| + <input type="checkbox"/> Footrest, long * (depth: 18 cm / 7.1"), tiltable and adjustable in angle | 9800008051010 |
| + <input type="checkbox"/> Footrest, short * (depth: 12 cm / 4.7"), tiltable and adjustable in angle | 9800008060010 |
| + <input type="checkbox"/> Seat depth extension * (+13 cm / 5.1") | 9800008082000 |
| + <input type="checkbox"/> Stand, front * (for swivel base) (recommended when using „Seatfix“, not necessary when using the footrest) | 9800008170010 |

Swivel base with adapter is required to use the + accessories.

- | | |
|--|----------------------|
| <input type="checkbox"/> Swivel base with adapter for accessories *** (incl. adapter plate) | 9800008202000 |
| <input type="checkbox"/> Back stabilisers (required when using swivel base) | 9800008161000 |
| + <input type="checkbox"/> Seatfix connection (for swivel base) ** option 1 (if vehicle is equipped with Isofix) | 9800008181000 |
| + <input type="checkbox"/> Stabilising bow, rear (for swivel base) option 2 (if the vehicle is NOT equipped with Isofix) | 9800008190000 |
| ↓ can be combined with: | |
| + <input type="checkbox"/> Footrest, long * (depth: 18 cm / 7.1"), tiltable and adjustable in angle | 9800008051010 |
| + <input type="checkbox"/> Footrest, short * (depth: 12 cm / 4.7"), tiltable and adjustable in angle | 9800008060010 |
| + <input type="checkbox"/> Stand, front * (for swivel base) (recommended when using „Seatfix“, not necessary when using the footrest) | 9800008170010 |

ACCESSORIES	item code
-------------	-----------

- | | |
|--|--------------------------------|
| <input type="checkbox"/> Abduction block | <i>included in basic model</i> |
| <input type="checkbox"/> Table, black | 9800008111012 |
| <input type="checkbox"/> Soft harness strap cover (pair) | 9800008151000 |
| <input type="checkbox"/> Chest clip | 9800004230000 |
| <hr/> | |
| <input type="checkbox"/> Seat wedge, inside | 9800008092000 |
| <input type="checkbox"/> Seat wedge, below ** (+10°) | 9800008104000 |
| <hr/> | |
| <input type="checkbox"/> Car seat protection mat incl. back | 9800008141000 |

Hint for more safety: With the anti-escape vest children remain belted securely in the car seat.

- | | |
|--|----------------------|
| <input type="checkbox"/> NapUp grey | 9800806550000 |
| <input type="checkbox"/> Anti-escape vest | 9800008260000 |

NOTE:

* Only footrest OR seat depth extension OR stand can be used.

** In some vehicles, the combination of seat wedge (below) and Seatfix connection cannot be used.

*** When using the swivel base, the Seatfix locking arms on the basic model may have to be removed depending on the vehicle type.