



**Rehab Safety Car Seat
for children**

**Restraint system
for special needs „S“**

ECE R 44/04 - Group I - III

**from 9 kg / 19.8 lb
up to 36 kg / 79.4 lb**

**Order e-mail:
info.export@thomashilfen.de**

**Status: 01.01.2024
valid until further notice**

Address

Customer no.:

Company:

Street:

Postal Code:

City:

Commission:

Contact person:

Phone:

E-mail:

Date:

BASIC MODEL	item code
-------------	-----------

- | | |
|--|----------------------|
| <input type="checkbox"/> RECARO Sport HERO Reha | 9800008802000 |
|--|----------------------|

Basic configuration: Seat and back with integrated 4-fold height adjustable headrest; integrated 5-point-belt with central adjuster; height-adjustable shoulder belt with HERO safety-system incl. crotch pad (small); advanced side protection (ASP); practical carrying handle; adjustable resting position; removable head, lateral trunk and side supports; removable covers, can be machine-washed at 30°C (86°F)


Measurements: seat depth: 31 cm / 12.2" (26 cm / 10.2" with supports) / seat with: front 31 cm / 12.2" (21 cm / 8.3" with supports), back 29 cm / 11.4" (20 cm / 7.9" with supports) / back height (seat): 55 - 73* cm / 21.7 - 28.7"* (* after removing guide rail) / shoulder height (child): 23 - 42 cm / 9.1 - 16.5"

ACCESSORIES	item code
-------------	-----------

Adapter or swivel base with adapter is necessary for the use of the + accessories.

- | | |
|---|----------------------|
| <input type="checkbox"/> Adapter for accessories (footrest, stand, stabilising bow, seatfix) ** | 9800008210000 |
| <input type="checkbox"/> Swivel base with adapter for accessories (footrest, stand, stabilising bow, seatfix) * / ** / *** | 9800008200000 |
| + <input type="checkbox"/> Footrest, short ** (depth: 12 cm / 4.7"), tiltable and adjustable in angle | 9800008061011 |
| + <input type="checkbox"/> Footrest pad (10 cm / 3.9") (only in combination with footrest, short) | 5005110501010 |
| <i>Recommended options for stability:</i> | |
| + <input type="checkbox"/> Seatfix connection (for swivel base) ** / **** <i>option 1</i> | 9800008181000 |
| + <input type="checkbox"/> Stabilising bow, rear (for swivel base) ** / **** <i>option 2</i> | 9800008190000 |
| + <input type="checkbox"/> Stand, front (for swivel base) ** (recommended with „seatfix connection“, not applicable when using footrest) | 9800008170010 |
| <hr/> | |
| <input type="checkbox"/> Seat wedge, below *** (+ 10°) | 9800008104000 |
| <input type="checkbox"/> Table | 9800008852000 |
| <hr/> | |
| <input type="checkbox"/> Belt extension (+ 20 cm / 7.9") incl. crotch belt (only for children in ECE group II + III, 18-36 kg) | 9800008240000 |
| <input type="checkbox"/> Crotch pad (large) | 9700910351000 |
| <input type="checkbox"/> Soft harness strap cover, grey (pair) (small pads as an alternative to standard belt pads) | 9800008151000 |
| <hr/> | |
| <input type="checkbox"/> Car seat protection mat incl. back | 9800008141000 |

Hint for more safety: With the anti-escape vest children remain belted securely in the car seat.

- | | | |
|--|--|----------------------|
| <input type="checkbox"/> NapUp grey |  | 9800806550000 |
| <input type="checkbox"/> Anti-escape vest |  | 9800008260000 |

NOTE:

* Swivel base cannot be used when fastening the RECARO Sport HERO Reha car seat with the 3-point-belt (only ECE group I: 9 kg / 19.8" up to 18 kg / 45.7") of the vehicle.
Recommended options for stability – You should choose a seatfix connection or a stabilising bow when ordering a swivel base!
 ** The adapter or swivel base with adapter must be ordered to use the footrest, stand, stabilising bow or seatfix connection.
 *** Swivel base and seat wedge (below) cannot be used in combination
 **** Seatfix connection and stabilising bow (rear) cannot be used in combination