

TECHNICAL MANUAL



EASyS Advantage Pediatric Wheelchair

- Models (seating system and chassis)
- Accessories for EASyS Advantage
- Accessories for A chassis Advantage and Q chassis

INDEX

EASyS Advantage ...

... is a pediatric wheelchair for children with safety seat.

You can find further information about Advantage on the Internet: www.thomashilfen.com or www.thomashilfen.us

Content	page
MODELS	EASyS Advantage seating system – size 1 and 2 3
	A chassis Advantage – size 1 and 2 4
	Q chassis. 4
MODELS standard features optional	Upholstery set 5
ACCESSORIES for EASyS Advantage seating system	Foot straps 6
	Calf strap 6
	Footrest-tube upholstery 6
	Abduction block 6
	Grip rail 7
	Grip rail, strong curved 7
	Lateral trunk support with chest belt 7
	Lateral trunk support 7
	Tray 8
	Accessory bag 8
	Frame padding 8
	2-point pelvic belt with front pull 9
	4-point pelvic belt with front pull 9
	Crotch strap 9
	Soft pelvic harness 9
	Butterfly/chest harness 10
	Butterflychest harness, flexible 10
	Seat waistcoat 10
	H-belt, padded 10
	Lateral head support 11
	Lateral head support, 45° 11
	Sun and rain cover with leg blanket 12
	Sun and rain cover 12
	Rain shield for sun and rain cover 12
	Transparent cover with leg blanket 13
	Sleeping bag (woven fur) 13
	Umbrella 13
ACCESSORIES for A chassis Advantage	Vent tray 14
	Pneumatic wheels 14
	Tip assist 14
	Transit tie down system 14
	Oxygen holder 14
ACCESSORIES for Q chassis	Adapter 15
	Push bar 15

SEATING SYSTEM



Advantage seating system	size 1	size 2
Item code	9800066098000	9800067098000
Seat depth	19 - 30 cm 7.5 - 11.8"	27 - 38 cm 10.6 - 15"
Seat width	19 - 31 cm 7.5 - 12.2"	25 - 37 cm 9.8 - 14.6"
Back height	41,5 - 68 cm 16.3 - 26.8"	55 - 74 cm 21.7 - 29.1"
Lower leg length	17.5 - 29 cm 6.9 - 11.4"	21 - 37 cm 8.3 - 14.6"
Hip angle = Back angle	90° - 180°	90° - 180°
Knee angle	90° - 180°	90° - 180°
Foot angle	78° - 100°	78° - 100°
Turning range of seat	180°	180°
Seat tilt-in-space in driving direction (on A chassis Advantage)	-10° to +35°	-10° to +35°
Seat tilt-in-space face-to-face (on A chassis Advantage)	-10° to +35°	-10° to +35°
Size of hip guide (w x h)	23 x 12 cm 9.1 x 4.7"	30 x 15 cm 11.8 x 5.9"
Size of footrest (w x d)	28 x 18 cm 11 x 7.1"	32 x 21 cm 12.6 x 8.3"
Weight (seat)	9 kg / 19.8 lb	10.4 kg / 22.9 lb
Max. user weight (seat)	35 kg / 77.2 lb	40 kg / 88.2 lb
Material	upholstery: Polyester – Trevira CS, 3-mesh- air-comfort-interlaced yarns (AGW)	
seat:	aluminium tubular frame, plastic assembly units	

Standard features

BSS (Balance Safety System) – safety seat; physiological pivots for the hip and knee joint; breathable seat and back upholstery as well as side padding (freely selectable) with reflective safety-seams; additional upholstery for more sitting comfort; fire retardant fabrics; lateral body guides; tiltable footrest, height and angle adjustable; brackets for sun and rain cover; continuously adjustable back height, hip angle, seat depth and seat width; easily accessible adjustments using star screws

CHASSIS



A chassis Advantage	size 1	size 2
Item code (swivel front wheels)	9800062014000	9800064014000
Item code (fixed front wheels)	9800062003000	9800064003000
Item code (Crash)	9800062024000	9800064024000
Item code (Model 22 Crash)	9800062025000	9800064025000
Seat tilt-in-space in driving direction	-10° to +35°	-10° to +35°
Seat tilt-in-space face-to-face	-10° to +35°	-10° to +35°
Turning range of seat	180°	180° / 18.1"
Seat height with seat unit	47 cm / 18.5"	46 cm / 18.1"
Push bar height	71.5 - 117 cm / 28.1 - 46.1"	71.5 - 117 cm / 28.1 - 46.1"
Turning circle	83 cm / 32.7"	99 cm / 32.7"
Overall dimensions (l x w x h) only A chassis Advantage	99.5 x 59.5 x 96 cm 39.2 x 23.4 x 37.8"	107 x 69 x 108 cm 42.1 x 27.2 x 42.5"
Folded dimensions (l x w x h) only A chassis Advantage	81 x 59.5 x 38 cm 31.9 x 23.4 x 15"	81 x 69 x 38 cm 31.9 x 27.2 x 15"
Overall dimensions (l x w x h) with seat unit EASyS Advantage	99.5 x 59.5 x 114 cm 39.2 x 23.4 x 44.9"	107 x 69 x 120 cm 42.1 x 27.2 x 47.2"
Folded dimensions (l x w x h) with seat unit EASyS Advantage	82 x 59.5 x 42.5 cm 32.3 x 23.4 x 16.7"	86 x 69 x 55 cm 33.9 x 27.2 x 21.7"
Wheel size swivel wheels (front)	18 cm / 7.1"	18 cm / 7.1"
Wheel size (rear)	25 cm / 9.8"	25 cm / 9.8"
Wheel size (Model 22 Crash front)	18.5 cm / 7.3"	18.5 cm / 7.3"
Wheel size (Model 22 Crash rear)	28.5 cm / 11.2"	28.5 cm / 11.2"
Weight (without seat unit)	9.8 kg / 21.6 lb	10.9 kg / 24 lb
Weight (without seat unit, Modell 22 Crash)	10.8 kg / 23.3 lb	11.9 kg / 26.2 lb
Max. load capacity basket	5 kg / 11 lb	5 kg / 11 lb
Max. load capacity (incl. seat unit)	50 kg / 110.2 lb	55 kg / 121.3 lb
Material	aluminium tubular frame, powder-coated, rim: PP (polypropylene); tires: composition rubber	

Standard features

Chassis: lightweight and stable aluminum chassis, foldable; adjustment of seat angle inclination as well as seat assembling and disassembling by central lever; included spring in rear frame; height adjustable push bar with biocompatible soft grip; foot brake; no-flat EVA-wheels; fixed front wheels or swivel front wheels incl. swivel lock; removable basket

The "Crash" model is also equipped with special attachment points for use as a seat in a motor vehicle equipped

Seat unit: see page 3

Q chassis	size 1	size 2
Item code chassis	9800066666000	9800066666000
Seat tilt-in-space	-7° to +31°	-7° to +31°
Seat height with seat unit	29 - 69 cm / 11.4 - 27.2"	29 - 69 cm / 11.4 - 27.2"
Adapter height	21 - 60 cm / 8.3 - 23.6"	21 - 60 cm / 8.3 - 23.6"
Size of swivel wheels (front)	100 mm / 3.9"	100 mm / 3.9"
Size of swivel wheels (rear with brake)	100 mm / 3.9"	100 mm / 3.9"
Total dimensions (only Q chassis) (l x w x h)	76 x 54 x 25 - 62 cm / 29.9 x 21.3 x 9.8 - 24.4"	76 x 54 x 25 - 62 cm / 29.9 x 21.3 x 9.8 - 24.4"
Weight (without seat unit)	13.5 kg / 29.8 lb	13.5 kg / 29.8 lb
Max. load capacity	50 kg / 110.2 lb	50 kg / 110.2 lb
Material	aluminium tubular frame, powder-coated, PU castors	


Standard features

Chassis: comfortable indoor chassis, easy height and seat tilting by gas pressure spring (foot pressure release), 4 swivel castors, 2 of them lockable (rear)


Seat unit: see page 3



Upholstery set	size 1	size 2
Item code	9800062470000	9800062490000
Color	red	red



Upholstery set	size 1	size 2
Item code	9800062480000	9800064480000
Color	blue	blue



Upholstery set	size 1	size 2
Item code	9800062490000	9800064490000
Color	green	green



Upholstery set	size 1	size 2
Item code	9800062460000	9800064460000
Color	magenta	magenta



Upholstery set	size 1	size 2
Item code	9800062450000	9800064450000
Color	gray	gray

ACCESSORIES



Foot straps	size 1	size 2
Item code	9800066482000	9800066482000
Belt length	87 cm / 34.3"	87 cm / 34.3"
Color	black	black



Calf strap	size 1	size 2
Item code	9800066355000	standard feature
Size (w x h)	29 x 10 x 1 cm 11.4 x 3.9 x 0.4"	–
Color	black	black



Footrest-tube upholstery	size 1	size 2
Item code	9800062300000	9800062300000
Size (w x h)	7 x 15 cm / 2.8 x 5.9"	7 x 15 cm / 2.8 x 5.9"
Color	gray	gray



Abduction block	size 1	size 2
Item code	9800066206000	9800067206000
Size (w x h x l)	6 x 9 x 10 cm 2.4 x 3.5 x 3.9"	7 x 10,5 x 11 cm 2.8 x 4.1 x 4.3"
Color	gray	gray



Grip rail	size 1	size 2
Item code	9800066215000	9800067215000
Height to seat surface	9 / 11,5 / 14 / 16,5 cm 3.5 / 4.5 / 5.5 / 6.5"	14 / 17 / 19 / 22 cm 5.5 / 6.7 / 7.5 / 8.7"
Width (inside dimension)	34 cm / 13.4"	38 cm / 13.4"
Color	gray	gray



Grip rail, strong curved	size 1
Item code	9800066215100
Height to seat surface	10.5 / 13 / 15,5 / 18.5 cm 4.1 / 5.1 / 6.1 / 7.3"
Width (inside dimension)	34 cm / 13.4"
Color	gray



Lateral trunk support, with chest belt	size 1	size 2
Item code	9800066332000	9800066332000
Supporting surface (h x d)	10 x 20 cm / 3.9 x 7.9"	10 x 20 cm / 3.9 x 7.9"
Color	gray	gray



Lateral trunk support	size 1	size 2
Item code	9800066346000	9800067346000
Size (h x d)	15 x 13 cm / 5.9 x 5.1"	17 x 15 cm / 6.7 x 5.9"
Color	gray	gray

ACCESSORIES



Tray	size 1	size 2
Item code	98000662810100	98000672810100
Size (w x d)	56 x 33.2 cm 5.9 x 8.7"	60 x 36.2 cm 23.6 x 14.3"
Height to seat surface	15 - 22 cm / 5.9 - 8.7"	15 - 22 cm / 5.9 - 8.7"
Edge height (inside)	2.5 cm / 1"	2.5 cm / 1"
Angle adjustment	approx. 27°	approx. 27°
Depth adjustment range	16 cm / 6.3"	19 cm / 7.5"
Cut-out (w x d)	41 x 14 cm 16.1 x 5.5"	45,5 x 17,5 cm 17.9 x 6.9"
Material thickness	8 mm / 0.3"	8 mm / 0.3"
Color tray plate	transparent / colorless	transparent / colorless



Accessory bag	size 1	size 2
Item code	9800066394100	9800066394100
Size (w x h x d)	30 x 24 x 7 cm 11.8 x 9.4 x 2.8"	30 x 24 x 7 cm 11.8 x 9.4 x 2.8"
Number of pockets	4 (2x small; 2x large)	4 (2x small; 2x large)
Pocket size small (w x h)	13 x 13 cm / 5.1 x 5.1"	13 x 13 cm / 5.1 x 5.1"
Pocket size large (w x h)	30 x 24 cm / 11.8 x 9.4"	30 x 24 cm / 11.8 x 9.4"
Color	black	black



Frame padding	size 1	size 2
Item code, black	9800066226000	9800066226000
Item code red-black	9800066227000	9800066227000
Outer dimensions (w x h)	(10 - 15) x 45 cm (3.9 - 5.9) x 17.7"	(10 - 15) x 45 cm (3.9 - 5.9) x 17.7"



	size 1	size 2
2-point pelvic belt with front pull		
Item code	9800084201000	9800084701000
Size of pad (b x h)	15 x 5 cm / 5.9 x 2"	17 x 5 cm / 6.7 x 2"
Belt length	130 cm / 51.2"	130 cm / 51.2"
Color	black	black



	size 1	size 2
4-point pelvic belt with front pull		
Item code	9800084211000	9800084711000
Size of pad (b x h)	15 x 8 cm / 5.9 x 3.1"	17 x 8 cm / 6.7 x 3.1"
Rear belt length	120 cm / 47.2"	120 cm / 47.2"
Front belt length	50 cm / 19.7"	50 cm / 19.7"
Color	black	black



	size 1	size 2
Crotch strap		
Item code	9800063411000	9800063411000
Belt length	approx. 18 - 33 cm approx. 7.1 - 13"	approx. 18 - 33 cm approx. 7.1 - 13"
Color	black	black
Note	can only be used in combination with pelvic harnesses	



	size 1	size 2
Soft pelvic harness		
Item code	9800066146000	9800067146000
Rear upholstery width	21 cm / 8.3"	26,5 cm / 8.3"
Belt length	45 cm / 17.7"	45 cm / 17.7"
Color	black	black



	size 1	size 2
5-point pelvic belt		
Item code	9800084800000	9800084800000
Size shoulder pads (w x l)	4.5 x 14.5 cm / 1.8 x 5.7"	4.5 x 14.5 cm / 1.8 x 5.7"
Shoulder belt length	86 cm / 33.9"	86 cm / 33.9"
Pelvis belt length	45 cm / 17.7"	45 cm / 17.7"
Color	black	black

ACCESSORIES



Butterfly/chest harness	size S	size M	size L
Item code	9800066166000	9800066176000	9800067176000
A	18.5 cm / 7.3"	25 cm / 9.8"	25.5 cm / 10"
B	30 cm / 11.8"	33.5 cm / 13.2"	45 cm / 17.7"
C	16 cm / 6.3"	18.5 cm / 7.3"	23 cm / 9.1"
D	14 cm / 5.5"	17 cm / 6.7"	22 cm / 8.7"
E	11.5 cm / 4.5"	10.5 cm / 4.1"	11.5 cm / 4.5"
F	18 cm / 7.1"	21.5 cm / 8.5"	25 cm / 9.8"



Butterfly/chest harness, flexible	size S	size M
Item code	9800063853000	9800065853000
A	12 cm / 4.7"	18 cm / 7.1"
B	21 cm / 8.3"	33 cm / 13"
C	11.5 cm / 4.5"	15 cm / 5.9"
D	9.5 cm / 3.7"	18 cm / 7.1"
E	6 cm / 2.4"	10 cm / 3.9"
F	21 cm / 8.3"	27 cm / 10.6"



Seat waistcoat	size S	size M	size L
Item code	9800066315100	9800066325100	9800067315100
Upholstered surface (w x l)	33 x 30 cm 13 x 11.8"	38 x 35 cm 15 x 13.8"	44 x 41 cm 17.3 x 16.1"
Belt length, top	65 cm / 25.6"	65 cm / 25.6"	65 cm / 25.6"
Color	black	black	black



H-belt, padded	size 1	size 2
Item code	9800066297000	9800067146000
Size Shoulder pads (w x l)	5 x 15 cm / 2 x 5.9"	5 x 15 cm / 2 x 5.9"
Shoulder belt length	60 cm / 23.6"	60 cm / 23.6"
Pelvis belt length	40 cm / 15.7"	40 cm / 15.7"
Color	black	black



Lateral head support	size 1	size 2
Item code	9800066246000	9800067246000
Size (h x d)	13 x 11 cm / 5.1 x 4.3"	15 x 14 cm / 5.9 x 5.5"
Color	gray	gray



Lateral head support, 45°	size 1	size 2
Item code	9800066231000	9800066231000
Size (h x d)	10 x 12 cm / 3.9 x 4.7"	10 x 12 cm / 3.9 x 4.7"
Color	gray	gray

ACCESSORIES



Sun and rain cover with leg blanket	size 1	size 2
Item code	9800066376000	9800066376000
Hood (w x h, inner dimension)	43 x 51 cm 16.9 x 20.1"	43 x 51 cm 16.9 x 20.1"
Hood pocket (w x h x d)	24 x 28 x 5 cm 9.4 x 11 x 2"	24 x 28 x 5 cm 9.4 x 11 x 2"
Leg blanket (w x l)	43 x 100 cm 16.9 x 39.4"	43 x 100 cm 16.9 x 39.4"
Color	black / light gray	black / light gray
Note	suitable up to max. 71 cm / 28" vertex height	



Sun and rain cover	size 1	size 2
Item code	98000663761000	9800066376100
Hood (w x h, inner dimension)	43 x 51 cm 16.9 x 20.1"	43 x 51 cm 16.9 x 20.1"
Hood pocket (w x h x d)	24 x 28 x 5 cm 9.4 x 11 x 2"	24 x 28 x 5 cm 9.4 x 11 x 2"
Color	black / light gray	black / light gray
Note	suitable up to max. 71 cm / 28" vertex height	



Rain shield	size 1	size 2
Item code	9800067390000	9800067390000
Hood (w x h, inner dimension)	43 x 51 cm 16.9 x 20.1"	43 x 51 cm 16.9 x 20.1"
Opening (w x h)	33 x 34 cm 13 x 13.4"	33 x 34 cm 13 x 13.4"
Leg blanket (w x l)	60 x 120 cm 23.6 x 47.2"	60 x 120 cm 23.6 x 47.2"
Color	transparent	transparent
Note	can only be used together with sun and rain cover	



Transparent cover with leg blanket	size 1	size 2
Item code	9800066381000	9800066381000
Hood (w x h, inner dimension)	43 x 51 cm 16.9 x 20.1"	43 x 51 cm 16.9 x 20.1"
Leg blanket (w x l)	43 x 100 cm 16.9 x 39.4"	43 x 100 cm 16.9 x 39.4"
Color	transparent	transparent
Note	suitable up to max. 71 cm / 28" vertex height	



Umbrella	size 1	size 2
Item code	9800065942000	9800065942000
Diameter Ø	65 cm / 25.6"	65 cm / 25.6"
Height	75 cm / 28.3"	75 cm / 28.3"
Color	grey	grey



Sleeping bag, woven fur	size 1	size 2
Item code	9800066185000	9800067185000
Bottom part of leg blanket (w x l)	35 x (68 - 80) cm 13.8 x (26.8 - 31.5)"	40 x (90 - 104) cm 15.7 x (35.4 - 40.9)"
Leg blanket (lower leg length)	31 cm / 12.2"	43 cm / 16.9"
Color	black / light gray	black / light gray

ACCESSORIES – A chassis Advantage



Vent tray	size 1	size 2
Item code	9800062941000	9800064941000
Size (w x h x l)	38 x 3 x 39 cm 15 x 1.2 x 15.4"	43 x 3 x 39 cm 16.9 x 1.2 x 15.4"

Pneumatic wheels	
Item code	9800066471100
Wheel size swivel wheels (front)	19.05 cm / 7.5"
Wheel size (rear)	25.4 cm / 10"

Tip assist	
Item code	9800067502000
Tread surface (w x l)	4 x 8 cm / 1.6 x 3.1"
Color	anthracite

Transit tie down system	
Item code	9700678391000
Note	For further information please contact us. +49 4761 886-63 or -68

Oxygen holder	
Item code	9800062970000
Tread area (w x l)	41,5 x 14,2 x 31 cm / 16.4 x 5.6 x 12.2"
Color	black

ACCESSORIES – Q chassis



Adapter	size 1	size 2
Item code	9800062951000	9800062951000
Note	always required when using Advantage seat unit	



Push bar	size 1	size 2
Item code	9800066601000	9800066601000
Height	52 - 100 cm (angle approx. 90°)	52 - 100 cm (angle approx. 90°)
Adjustment angle	0 - 110°	0 - 110°
Color	black	black



Member of the International
Support Association for the
Rehabilitation of Children and Youths e.V.



For customers worldwide (except U.S.):

Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG
Walkmühlenstraße 1 · 27432 Bremervörde · Germany
Phone: +49 (0) 4761 886-63 or -68 · Fax: +49 (0) 4761 886-19
info@thomashilfen.de · www.thomashilfen.com

For customers in the U.S.:

Thomashilfen North America
14221 29th St E, #101– Sumner, WA 98390
Phone: 206 763 0754 – Toll free: 866.870.2122
info@thomashilfen.de · www.thomashilfen.us