



NEW

EASyS Advantage

Maximum driving comfort
for an easy-going day


Thomashilfen

More comfort, more safety: This is **NEW** ...

” Sitting for long periods can be exhausting, but sometimes it is necessary. You will be thankful when your rehab stroller is perfectly adjustable with cozy, comfortable seating. “



Easy seat assembly and rapid change of viewing direction – without tools:

- forward facing
- rear facing



180° lying position for children who need it



A single release lever on the back

- seat tilt-in-space
- separates the seat from the base frame



Transportation lock to prevent unwanted unfolding, e.g. when lifting the stroller into the car

Extremely thick cushions for great comfort.

No tools needed for seat width and depth adjustments

Choose between 5 upholstery styles:

				
red/ grey	blue/ grey	green/ grey	magenta/ grey	grey/ grey

Cushions are breathable and washable

Chassis colour: black 



New sportier footrest



Lever for easy adjustment of the knee angle at the pivot point

Much more storage space available under the seat:

- washable basket (inclusive)
- vent tray (optional)



Soft suspension absorbs vibrations and ensures a good riding pleasure without compromising steering



Convenient brake, clearly visible in red

EASyS Advantage

Maximum driving comfort for an easy-going day

Accessories (selection)



Handy accessory plate for medical devices, e.g. oxygen tank or suction equipment



Grip rail and flexible butterfly/chest harness for a comfortable and secure position of the upper body



Abduction block for guiding the thigh



Flexible lateral trunk supports and a 4-point pelvis belt for positioning the upper body and the pelvis



Comfortable sleeping bag, padded inside or with woven fur for cold days



Handy tip assist to facilitate curbs and other obstacles



Accessory bag on the back for all the little necessities

You will find the total overview of all rehab accessories in the order form.

Technical data (seat)	EASyS Advantage 1	Advantage 2
Seat depth	19 - 30 cm 7.5 - 11.8"	27 - 38 cm 10.6 - 15"
Seat width	19 - 31 cm 7.5 - 12.2"	25 - 37 cm 9.8 - 14.6"
Seat height	47 cm / 18.5"	46 cm / 18.1"
Back height	41,5 - 68 cm 16.3 - 26.7"	55 - 74 cm 21.7 - 29.1"
Lower leg length	17,5 - 29 cm 6.9 - 11.4"	21 - 37 cm 8.3 - 14.6"
Hip angle	90° - 180°	90° - 180°
Knee angle	90° - 180°	90° - 180°
Foot angle	78° - 100°	78° - 100°
Turning range of seat	180°	180°
Size of hip guide (w x h)	23 x 12 cm 9.1 x 4.7"	30 x 15 cm 11.8 - 5.9"
Size of foot rest (w x d)	28 x 18 cm 11 x 7.1"	32 x 21 cm 12.6 - 8.3"
Weight seat unit	9 kg / 19.8 lbs	10.4 kg / 22.9 lbs
max. user weight seat unit	35 kg / 77.2 lbs	40 kg / 88.2 lbs

Technical data	A1 chassis	A2 chassis
Seat tilt-in-space driving direction	-10° to +35°	-10° to +35°
Seat tilt-in-space face-to-face	-10° to +35°	-10° to +35°
Push bar height	71.5 - 117 cm 28.2 - 46.1"	71.5 - 117 cm 28.2 - 46.1"
Wheel size swivel wheels (front)	18 cm / 7.1"	18 cm / 7.1"
Wheel size (rear)	25 cm / 10"	25 cm / 10"
Weight chassis	9.8 kg / 21.6 lbs	10.9 kg / 24 lbs
max. carrying capacity basket	5 kg / 11 lbs	5 kg / 11 lbs
max. carrying capacity (incl. seat)	50 kg / 110.2 lbs	55 kg / 121.3 lbs

Technical data stroller	EASyS Advantage 1	Advantage 2
Overall dimensions (l x w x h)	99.5 x 59.5 x 114 cm 39.2 x 23.4 x 48.9"	107 x 69 x 120 cm 42.1 x 27.2 x 47.2"
Folded dimensions (l x w x h)	82 x 59.5 x 42.5 cm 32.3 x 23.4 x 16.7"	86 x 69 x 55 cm 33.9 x 27.2 x 21.7"