



**Rehab stroller for children**

**Order e-mail:**  
info.export@thomashilfen.de

**Status: 15.02.2024**  
valid until further notice

**Address** **Customer no.:**

Company: \_\_\_\_\_

Street: \_\_\_\_\_

Postal Code: \_\_\_\_\_

City: \_\_\_\_\_

Commission: \_\_\_\_\_

Contact person: \_\_\_\_\_

Phone: \_\_\_\_\_

E-mail: \_\_\_\_\_

Date: \_\_\_\_\_

BASIC MODEL	item code size 1	item code size 2
-------------	------------------	------------------



**Swiftly stroller**

**9800068065000**

**Standard features Swiftly:** lightweight and foldable chassis with integrated seat unit; crash eyelets; 5-point pelvic harness, 25 mm; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels with wheel lock; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover

**Measurements:**

- seat depth: 22 cm / 8.7" (17 cm / 6.9" with seat minimizer) - 28,5 cm / 11.2"
- seat width: 34 cm / 13.4" (23 cm / 9.1" with seat minimizer)
- back height: 62 cm / 24.4"
- min. shoulder height: 37 cm / 14.6" (26 cm / 10.2" with seat minimizer)



**Swiftly 2 stroller**

**9800068080010**

**Standard features Swiftly 2:** lightweight and foldable chassis with integrated seat unit; crash eyelets; 5-point pelvic harness, 25 mm; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels with wheel lock; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover; folding lock

**Measurements:**

- seat depth: 34 cm / 13.4" (24.5 cm / 9.6" with seat minimizer) - 40 cm / 15.7"
- seat width: 39.5 cm / 15.6" (29.5 cm / 11.6" with seat minimizer)
- back height: 70 cm / 27.6"
- min. shoulder height: 47 cm / 18.5" (30 cm / 11.8" with seat minimizer)

 accessories see page 2

# Swifty Order form (accessories)

Customer no.:

Company:

Contact person:

**Status: 15.02.2024**  
valid until further notice

**Order e-mail:**  
info.export@thomashilfen.de

ACCESSORIES	item code size 1	item code size 2
<b>FEET / LEGS</b>		
Foot straps black (pair)	<input type="checkbox"/> 9800066482000	<input type="checkbox"/> 9800066482000
Basket	<input type="checkbox"/> 9800068131000	<input type="checkbox"/> 9800068350000
Tip assist	<input type="checkbox"/> 9800068291000	<input type="checkbox"/> 9800068291000
<b>SEAT</b>		
Abduction block, grey (not possible in combination with crotch strap)	<input type="checkbox"/> 9800066206000	<input type="checkbox"/> 9800067206000
Grip rail, grey	<input type="checkbox"/> 9800068331000	<input type="checkbox"/> 9800067215000
Tray, transparent	<input type="checkbox"/> 9800068220000	<input type="checkbox"/> 9800067281010
Lateral trunk support, grey	<input type="checkbox"/> 9800066346000	<input type="checkbox"/> 9800067346000
Lateral trunk support with chest belt	<input type="checkbox"/> 9800066332000	<input type="checkbox"/> 9800066332000
Lateral seat upholstery, red	<input type="checkbox"/> 5005100228000	<input type="checkbox"/> 9700442011000
Frame padding black	<input type="checkbox"/> 9800066226000	<input type="checkbox"/> 9800066226000
Frame padding, red-black	<input type="checkbox"/> 9800066227000	<input type="checkbox"/> 9800066227000
<b>SEAT MINIMIZER</b>		
Seat minimizer, complete set, black	<input type="checkbox"/> 9800068340000	<input type="checkbox"/> 9800068360000
<i>The following components of the set can also be ordered separately:</i>		
Head pillow, black (Seat minimizer)	<input type="checkbox"/> 5005100209005	-
Backrest upholstery, black (Seat minimizer)	<input type="checkbox"/> 5005100219005	-
Pelvic upholstery, black (Seat minimizer)	<input type="checkbox"/> 5005100229005	-
Seat width reducer, black (Seat minimizer) (seat width: 29.5 cm / 11.6")	-	<input type="checkbox"/> 9800068410000
Footpad (Seat minimizer)	<input type="checkbox"/> 5005100239005	-
Seat minimizer, complete set, red	<input type="checkbox"/> 9800068341000	-
<i>The following components of the set can also be ordered separately:</i>		
Head pillow, red (Seat minimizer)	<input type="checkbox"/> 5005100225000	<input type="checkbox"/> 9800068420000
Backrest upholstery, red (Seat minimizer)	<input type="checkbox"/> 5005100226000	-
Pelvic upholstery, red (Seat minimizer)	<input type="checkbox"/> 5005100227000	-
Seat width reducer, red (Seat minimizer) (seat width: 29.5 cm / 11.6")	-	<input type="checkbox"/> 9800068411000
Footpad (Seat minimizer)	<input type="checkbox"/> 5005100239005	-
<b>HARNESSES</b>		
2-point pelvic harness with front-pull	<input type="checkbox"/> 9800084201000	<input type="checkbox"/> 9800084701000
4-point pelvic harness with front-pull	-	<input type="checkbox"/> 9800084711000
Crotch strap (fits for 2-point pelvic belt, 4-point pelvic belt)	<input type="checkbox"/> 9800063411000	<input type="checkbox"/> 9800063411000
Soft pelvic harness, black	<input type="checkbox"/> 9800066146000	<input type="checkbox"/> 9800067146000
Butterfly/chest harness, black (S)	<input type="checkbox"/> 9800066166000	-
Butterfly/chest harness, black (M)	<input type="checkbox"/> 9800066176000	<input type="checkbox"/> 9800066176000
Butterfly/chest harness, black (L)	-	<input type="checkbox"/> 9800067176000
Butterfly/chest harness, flexible (S)	<input type="checkbox"/> 9800063853000	-
Butterfly/chest harness, flexible (M)	<input type="checkbox"/> 9800065853000	<input type="checkbox"/> 9800065853000
Seat waistcoat, black (S)	<input type="checkbox"/> 9800066315100	-
Seat waistcoat, black (M)	<input type="checkbox"/> 9800066325100	<input type="checkbox"/> 9800066325100
Seat waistcoat, black (L)	-	<input type="checkbox"/> 9800067315100
<b>HEAD</b>		
Head pillow, black	<input type="checkbox"/> 9800068240000	<input type="checkbox"/> 9800068370000
Head pillow, red	-	<input type="checkbox"/> 9800068371000
<b>WEATHER PROTECTION</b>		
Sleeping bag, woven fur (body contour), winter	<input type="checkbox"/> 9800066185000	<input type="checkbox"/> 9800067185000
Rain cover, transparent, with leg blanket	<input type="checkbox"/> 9800066381000	<input type="checkbox"/> 9800066381000
Sun and rain cover with leg blanket	<input type="checkbox"/> 9800066376000	<input type="checkbox"/> 9800066376000
Sun and rain cover without leg blanket	<input type="checkbox"/> 9800066376100	<input type="checkbox"/> 9800066376100
Rain shield for sun and rain cover, transparent (can be used only in combination with „sun and rain cover“)	<input type="checkbox"/> 9800067390000	<input type="checkbox"/> 9800067391000
Umbrella with UV protection	<input type="checkbox"/> 9800065942000	<input type="checkbox"/> 9800065942000
<b>CHASSIS</b>		
Pneumatic wheels (upgrade kit)	<input type="checkbox"/> 9800066470100	<input type="checkbox"/> 9800066711000
Oxygen holder	<input type="checkbox"/> 9800062970000	<input type="checkbox"/> 9800062970000