

Recommended for
walking therapy!

trio

■ walker ■ ■



Now with
adjustable width!


Thomashilfen

3 advantages of the walker

- + puncture safe „Flexel“-tyres
- + high stability and tilt safety
- + tested according to DIN EN ISO 11199-3-2005

Further advantages:

- high stability and tilt safety, tested according to DIN EN ISO 11199-3-2005
- comfortable arm rest, height and width adjustable.
- adjustable without tools
- ideal for assembly of purpose built items (e.g. pelottes)
- height adjustable armpit supports (optional)
- single-hand brake possible (accessory)
- lockable hand brakes
- small turn radius
- walking practise inside



Height adjustable and tilting armpit supports (optional).



Comfortable arm rest, height and width adjustable. (trio B)



Single-hand brake retrofit kit (optional).

Product description	Item code	length	width	height (arm rest)	wheel-Ø	max. weight capacity
Walker trio basic	9800020501100	90 cm 35.4"	71 - 77 cm 28 - 30.3"	90 - 130 cm 35.4 - 51.2"	200 x 50 mm 7.9 x 2"	150 kg 331 lb.
Walker trio B (+ width adjustment)	9800020451100	90 cm 35.4"	71 - 77 cm 28 - 30.3"	90 - 130 cm 35.4 - 51.2"	200 x 50 mm 7.9 x 2"	150 kg 331 lb.
Single-hand brake retrofit kit (trio)	9800018260003	-	-	-	-	-
Height adjustable and tilting armpit supports	9700182200000	-	-	-	-	-