



**Therapy Chair for children
for active, unrestraint sitting**

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:







Date:

Contact:

**Order fax:
+49 (0)4761 886-19**

Status: 01.01.2015
valid until further notice

BASIC MODEL	item code
-------------	-----------

- | | | | |
|--------------------------|---|---|-------------|
| <input type="checkbox"/> |  | <p>SiiS M – for the smallest
perfectly suited for the use in kindergartens
 seat height: 32 – 40 cm (12.6 – 15.8")*
 Features: • adjustable impulse generators and foot brake as standard
 • small seat block / small iliac crest support</p> | 8301 |
| <input type="checkbox"/> |  | <p>SiiS 1
narrower base frame
 seat width: 19 cm (7.5"); seat height: 42 – 55 cm (16.5 – 21.7")*;
 total dimensions (w x h x d): 54 x 88 x 54 cm (21.3 x 34.7 x 21.3")</p> | 8300 |
| <input type="checkbox"/> |  | <p>SiiS 2
wider base frame for more legroom
 seat width: 19 cm (7.5"); seat height: 42 – 55 cm (16.5 – 21.7")*;
 total dimensions (w x h x d): 60 x 88 x 56 cm (23.6 x 34.7 x 22")</p> | 8305 |
| <input type="checkbox"/> |  | <p>SiiS H – particularly suitable for persons with multiple sclerosis
perfectly suited for the use in sheltered workshops
 seat height: 53 – 72 cm (20.9 – 31.5")*
 Features: • less abduction
 • adjustable impulse generators as standard
 • longer brake lever</p> | 8307 |
| <input type="checkbox"/> |  | <p>SiiS S
 seat width: 16 cm (6.3"); seat height: 42 – 55 cm (16.5 – 21.7")*
 Features: • narrower seat surface</p> | 8308 |
| <input type="checkbox"/> |  | <p>SiiS P
perfectly suited for the use in sheltered workshops
 seat width: 22 cm (8.7"); seat height: 49 – 68 cm (19.3 – 26.8")*
 Features: • wider seat surface
 • higher seat height</p> | 8306 |

* 5 cm / 2" thereof by moving the gas spring unit upwards

Standard features:

underframe with four swivelling castors; ergonomic seat block; iliac crest support; 2 impulse generators; push bar (not standard at SiiS M); removable lever brake (not standard at SiiS M)

ACCESSORIES	item code
-------------	-----------

- | | | |
|--------------------------|--|---------------|
| <input type="checkbox"/> | <p>front wheel, lockable (piece)</p> | 8321 |
| <input type="checkbox"/> | <p>footbrake, right (standard for SiiS M)</p> | E61985 |
| <input type="checkbox"/> | <p>impulse generators, depth adjustable (piece) (standard for SiiS M and H)</p> | 8320 |