



**Therapy chair
with mobility seat
for children**

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

Order fax:
+49 (0)4761 886-19

Status: 10.02.2012
valid until further notice

BASIC MODEL	item code size 1	item code size 2
-------------	------------------	------------------

ThevoTherapy

8400

8450

Standard features:

ThevoTherapy seating system:

- flexible wing elements and wing suspensions in the back
- stable seat with anatomic guidance
- physiological pivots for the hip- and knee joint
- gas spring operated hip angle adjustment
- breathable seat and back upholstery
- fire retardant fabrics
- femoral support
- height and angle adjustable footrest

chassis:

- comfortable and solid indoor base
- height and seat tilt-in-space adjustment with gas spring
- easy foot engaged adjustment
- 4 swivelling castors for more comfort, 2 of them lockable

ACCESSORIES	item code size 1	item code size 2
-------------	------------------	------------------

CHASSIS

Push bar (cannot be used in combination with push handle)

6660

6660

FEET

Foot straps, black

6648/1

6648/1

Calf strap, black

6635/5

standard

SEAT

Abduction block, black

6620/5

6720/5

Tray, fold-away

8425

8475

Side plate, high

8423

8473

Pelvic support

8432

8479

Armrests, pair

8431

8478

Impulse sensors

8422

8472

Moldable aluminum lath

6384

6584

Upholstery for impulse sensors (must be ordered for moldable aluminium lath)

500512047

500512048

Lateral trunk support plane

8426

8476

Lateral trunk support with chest belt

8427

8427

Push handle (cannot be used in combination with push bar at the Therapy chassis)

8439

8439

HARNESSES

2-point pelvic belt with front pull

8420

8470

4-point pelvic belt with front pull

8421

8471

Soft pelvic harness, black

6614/51

6714/51

Butterfly / chest harness, flexible (S)

6385/2

-

Butterfly / chest harness, flexible (M)

6585/2

6585/2

Belt bow

6371/2

6571/2

Note: The belt bow has to be ordered in combination with the butterfly/chest harness, otherwise the mounting of the belts is not possible.

HEAD

Headrest, neck support

8429

8429

Headrest, shell

8424

8424

Headrest height reducer

6569

6569