



**Rehab Safety Car Seat  
for children, foldable**

**Restraint system  
for special needs „S“**

**ECE R 44/04 - Group II - III**

**from 15 kg / 33.1 lbs  
to 36 kg / 79.4 lbs**

**Order fax:  
+49 (0)4761 886-19**

**Status: 17.04.2019**  
valid until further notice

**Address**

Company: \_\_\_\_\_

Street: \_\_\_\_\_






ZIP Code / Place: \_\_\_\_\_

Phone: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

BASIC MODEL	item code
-------------	-----------

- |                          |   |  |            |
|--------------------------|---|--|------------|
| <input type="checkbox"/> |    | Navigator Car Seat, grey <span style="display: inline-block; width: 20px; height: 10px; background-color: #555; border: 1px solid #000; margin-left: 5px;"></span>       | <b>600</b> |
| <input type="checkbox"/> |   | Navigator Car Seat, red <span style="display: inline-block; width: 20px; height: 10px; background-color: #f00; border: 1px solid #000; margin-left: 5px;"></span>        | <b>601</b> |
| <input type="checkbox"/> |  | Navigator Car Seat, blue <span style="display: inline-block; width: 20px; height: 10px; background-color: #00a0e3; border: 1px solid #000; margin-left: 5px;"></span>    | <b>602</b> |
| <input type="checkbox"/> |  | Navigator Car Seat, green <span style="display: inline-block; width: 20px; height: 10px; background-color: #90ee90; border: 1px solid #000; margin-left: 5px;"></span>   | <b>603</b> |
| <input type="checkbox"/> |  | Navigator Car Seat, magenta <span style="display: inline-block; width: 20px; height: 10px; background-color: #e03399; border: 1px solid #000; margin-left: 5px;"></span> | <b>604</b> |

**Basic configuration:**

Seat and back with 7-fold height adjustable headrest; complete seat foldable; seat cushion with abduction block (2-fold adjustable); back cushion with side pads; head cushion; additional head cushions (left, right); 5-point positioning belt; chest clip; cupholder (applicable right or left); fastening belt for use outside the car; removable covers, can be machine-washed at 40°C (104°F)

**Measurements:** seat depth: 32 cm / 12.6" // seat width: front: 26 cm / 10.2", back: 26 cm / 10.2" // back height (seat): 61 – 71 cm / 24 – 28"

ACCESSORIES	item code
-------------	-----------

**Hint for more safety:** With the anti-escape vest children remain belted securely in the car seat.

- |                          |   |            |
|--------------------------|---|------------|
| <input type="checkbox"/> | Anti-escape vest, size 1 (for children with shoulder height <b>under</b> 31 cm / 12.2") | <b>822</b> |
| <input type="checkbox"/> | Anti-escape vest, size 2 (for children with shoulder height <b>over</b> 31 cm / 12.2")  | <b>823</b> |
| <input type="checkbox"/> | Transportation bag  | <b>605</b> |
| <input type="checkbox"/> | Car seat protection mat   | <b>813</b> |
| <input type="checkbox"/> | Car seat protection mat incl. back  | <b>814</b> |

