



Positioning harness for children and adults with lockable buckles for use in vehicles

Body weight: 30 - 80 kg / 66.1 - 176.4 lb

Body height: 130 - 185 cm / 51.2 - 72.8"

(up to 150 cm / 59.1" always use with seat shell / booster seat)

Order fax: +49 (0)4761 886-19

Status: 11.08.2020
valid until further notice

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

BASIC MODEL	item code
-------------	-----------

- | | |
|--|-----|
| <input type="checkbox"/> Octopus positioning harness | 690 |
|--|-----|

Basic configuration:

Base back part; 7-point seat vest with belt buckles; shoulder straps; pelvis belts; crotch strap; 7x lockable side release buckles; keys for belt buckles; emergency hammer with seatbelt cutter; water-repellent cushions

Measurements cushion (W x H): 7-point seat vest: 34 x 42 cm / 13.4 x 16.5" // Base back part: 34 x 84 cm / 13.4 x 33.1"

ACCESSORIES	item code
-------------	-----------


- | | | |
|--------------------------|--|-----|
| <input type="checkbox"/> |  <p>Belt extension</p> <p><i>(required for vehicle seats with firmly connected seat and backrest or fixed backrest, e.g. in buses)</i></p> | 691 |
|--------------------------|--|-----|

- | | | |
|--------------------------|---|-----|
| <input type="checkbox"/> |  <p>Base back part for positioning harness*</p> <p><i>(Features: Base back part; shoulder straps; pelvis belts; crotch strap; keys for belt buckles; emergency hammer with seatbelt cutter; water-repellent padding)</i></p> | 692 |
|--------------------------|---|-----|

- | | | |
|--------------------------|--|--|
| <input type="checkbox"/> |  <p>Headrest Cardiff</p> |  <p>80627</p> |
|--------------------------|--|--|

TIP:

* Recommended for a quick change of vehicle, e.g. school taxi, transport for the disabled, second car. The base back part for the positioning harness remains fully assembled in the vehicle, the 7-point seat vest changes vehicle quickly together with the user.

 **Please note:** The Octopus positioning harness is not classified as a safety harness, but as a medical device for positioning. To comply with safety regulations, the vehicle safety harness **MUST ALWAYS** be worn over the Octopus positioning harness. This combination provides additional safety for the user.