



**Rehab Safety Car Seat  
for children**

**Restraint system  
for special needs „S“**

**ECE R 44/04 - Group II - III**

**from 15 kg / 33.1 lbs  
tested up to 50 kg / 110.2 lbs**

**Order fax:  
+49 (0)4761 886-19**

**Status: 29.03.2017**  
valid until further notice

**Address**

Company: \_\_\_\_\_

Street: \_\_\_\_\_

ZIP Code / Place: \_\_\_\_\_

Phone: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

**Customer no.:** \_\_\_\_\_

BASIC MODEL	item code
-------------	-----------

- |  |            |
|--|------------|
| <input type="checkbox"/> <b>RECARO Monza Nova 2 Reha</b> | <b>800</b> |
|--|------------|

Seat and back with integrated 11-fold height adjustable headrest; stepless neck support of the headrest by an integrated air cushion; side impact protection by energy absorbing elements; removable head, lateral trunk and side supports; fixed armrests; 5-point positioning belt; crotch pad; removable abduction block; integrated speakers in the headrest for mp3-player use, integrated mp3-player pocket; RECARO Air Circulation System for optimum air circulation; ventilated foam; removable covers, can be machine-washed at 30°C (86°F)

- |  |            |
|--|------------|
| <input type="checkbox"/> <b>RECARO Monza Nova 2 Seatfix Reha</b> | <b>801</b> |
|--|------------|

Same features as RECARO Monza Nova 2 Reha; including integrated Isofix arms which can be fixed to the location points of vehicles which are equipped with Isofix

**NOTE:** Model Seatfix and swivel base cannot be used in combination!

**Measurements:** seat depth: 30 cm / 11.8" (43 cm / 16.9" with seat depth extension) / seat width: front 29 cm / 11.4", back 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside) / back height: 60 - 77 cm / 23.6 - 30.3"

**Hint for more safety:** With the anti-escape vest children remain belted securely in the car seat.

- |  |            |
|--|------------|
| <input type="checkbox"/> <b>Anti-escape vest, size 2</b> (for children with shoulder height <b>over</b> 31 cm / 12.2") | <b>823</b> |
|--|------------|

ACCESSORIES	item code
-------------	-----------

- |  |                                |
|--|--------------------------------|
| <input checked="" type="checkbox"/> <b>Abduction block</b>   | <i>included in basic model</i> |
| <input type="checkbox"/> <b>Footrest, long*</b> (depth: 18 cm / 7.1"), tiltable and adjustable in angle                                  | <b>805</b>                     |
| <input type="checkbox"/> <b>Footrest, short*</b> (depth: 12 cm / 4.7"), tiltable and adjustable in angle                                 | <b>806</b>                     |
| <input type="checkbox"/> <b>Swivel base with footrest adapter**</b> (incl. adapter plate)  | <b>820/1</b>                   |
| <input type="checkbox"/> <b>Back stabilisers</b> (required when using swivel base)   | <b>816</b>                     |
| <b>Recommended options for stability:</b>  |                                |
| <input type="checkbox"/> <b>Seatfix connection</b> (for swivel base) <i>option 1</i>   | <b>818</b>                     |
| <input type="checkbox"/> <b>Stabilising bow, rear</b> (for swivel base) <i>option 2</i>  | <b>819</b>                     |
| <input type="checkbox"/> <b>Stand, front</b> (for swivel base) (recommended when using „Seatfix“, not necessary when using the footrest) | <b>817</b>                     |

- |   |            |
|---|------------|
| <input type="checkbox"/> <b>Seat depth extension*</b> (+13 cm / 5.1") | <b>808</b> |
| <input type="checkbox"/> <b>Seat wedge, inside</b>                    | <b>809</b> |
| <input type="checkbox"/> <b>Seat wedge, below**</b> (+10°)            | <b>810</b> |
| <input type="checkbox"/> <b>Table</b>                                 | <b>811</b> |

- |  |            |
|--|------------|
| <input type="checkbox"/> <b>Soft harness strap cover</b> (pair)    | <b>815</b> |
| <input type="checkbox"/> <b>Car seat protection mat</b>            | <b>813</b> |
| <input type="checkbox"/> <b>Car seat protection mat incl. back</b> | <b>814</b> |

**NOTE:**

**Recommended options for stability – You should choose 819 or 818 when ordering a swivel base 820/1!**

\* Footrest and seat depth extension cannot be used in combination. Exception: The swivel base is also used.

\*\* Swivel base and seat wedge, below cannot be used in combination.