



Rehab Safety Car Seat
for children

Restraint system
for special needs „S“

ECE R 44/04 - Group I - III

from 9 kg / 19.8 lb
up to 36 kg / 79.4 lb

Order fax:
+49 (0)4761 886-19

Status: 16.11.2018
valid until further notice

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

BASIC MODEL	item code
-------------	-----------


- | | |
|--|------------|
| <input type="checkbox"/> RECARO Sport HERO Reha | 880 |
|--|------------|

Basic configuration: Seat and back with integrated 4-fold height adjustable headrest; integrated 5-point-belt with central adjuster; height-adjustable shoulder belt with HERO safety-system incl. crotch pad (small); advanced side protection (ASP); practical carrying handle; adjustable resting position; removable head, lateral trunk and side supports; removable covers, can be machine-washed at 30°C (86°F)

Measurements: seat depth: 31 cm / 12.2" (26 cm / 10.2" with supports) / seat with: front 31 cm / 12.2" (21 cm / 8.3" with supports), back 29 cm / 11.4" (20 cm / 7.9" with supports) / back height (seat): 55 - 73* cm / 21.7 - 28.7"* (* after removing guide rail) / shoulder height (child): 23 - 42 cm / 9.1 - 16.5"

ACCESSORIES	item code
-------------	-----------

Hint for more safety: With the anti-escape vest children remain belted securely in the car seat.

- | | |
|---|------------|
| <input type="checkbox"/> Anti-escape vest, size 1 (for children with shoulder height under 31 cm / 12.2") | 822 |
| <input type="checkbox"/> Anti-escape vest, size 2 (for children with shoulder height over 31 cm / 12.2") | 823 |
- 
- | | |
|---|---------------|
| <input type="checkbox"/> Footrest, short** (depth: 12 cm / 4.7"), tiltable and adjustable in angle | 806 |
| <input type="checkbox"/> Footrest pad (10 cm / 3.9") | E11050 |
| <input type="checkbox"/> Footrest adapter** | 821 |
| <input type="checkbox"/> Turning plate with footrest adapter*** | 820 |
| <i>Recommended options for stability:</i> | |
| <input type="checkbox"/> Seatfix connection (for turning plate)**** <i>option 1</i> | 818 |
| <input type="checkbox"/> Stabilising bow, rear (for turning plate)**** <i>option 2</i> | 819 |
| <input type="checkbox"/> Stand, front (for turning plate) | 817 |

- | | |
|--|------------|
| <input type="checkbox"/> Seat wedge, below*** (+ 10°) | 810 |
| <input type="checkbox"/> Table | 885 |

- | | |
|---|---------------|
| <input type="checkbox"/> Belt extension (+ 20 cm / 7.9") incl. crotch belt (only for children in ECE group II + III, 18-36 kg) | 824 |
| <input type="checkbox"/> Crotch pad (large) | E91035 |

- | | |
|--|------------|
| <input type="checkbox"/> Car seat protection mat | 813 |
| <input type="checkbox"/> Car seat protection mat incl. back | 814 |

NOTE:

Turning plate cannot be used when fastening the RECARO Sport HERO Reha car seat with the 3-point-belt (only ECE group I: 9 kg / 19.8" up to 18 kg / 45.7") of the vehicle.

Recommended options for stability – You should choose 819 or 818 when ordering a turning plate 820!

** The adapter must be ordered to use the footrest

*** Turning plate and seat wedge (below) cannot be used in combination

**** Seatfix connection and stabilising bow (rear) cannot be used in combination