

RECARO – the car seat that meets the highest safety standards; particularly growing with the child



9 - 36 kg
19.8 - 79.4 lb
 body weight



Measurements & weights RECARO Sport HERO Reha

Seat depth	31 cm / 12.2" (26 cm / 10.2" with seat pads)
Seat width	front 31 cm / 12.2" (21 cm / 8.3" with seat pads, 27 cm / 10.6" with seat minimizer) rear 29 cm / 11.4" (20 cm / 7.9" with seat pads, 25 cm / 9.8" with seat minimizer)
Back height (seat)	53 - 73 cm / 20.9 - 28.7"
Shoulder height (child)	23 - 42 cm / 9.1 - 16.5"
Overall dimensions (w x h x d)	54 x 55 (max. 81) x 52 cm / 21.3 x 21.7 (max. 31.9) x 20.5"
Weight of seat	approx. 8 kg / 17.7 lb
Max. carrying capacity	9 - 36 kg / 19.8 - 79.4 lb


RECARO & **Thomashilfen**
 A strong partnership

Car seat **RECARO Sport HERO Reha**

Car seat **RECARO Sport HERO Reha**

- Easy placing and removing of the children thanks to the turning plate (for ECE group II – III)
- Integrated seat tilt-in-space to rest and sleep position
- Pad system as standard
- Advanced Side Protection: Lateral energy absorbers reduce forces in case of a side impact
- 5-point safety belt with practical central adjuster (for ECE group I) incl. **innovative HERO Safety System**
- 5-point positioning belt optional (for ECE group II – III)
- ECE group I – III, approved acc. to current test standard ECE R 44/04

Functions



5-point harness with HERO technology



Fourfold height adjustable headrest



Easy and quick setting to rest position (for children of ECE group I only)



Practical carry handle

Accessories (selection)




Turning plate with Seatfix connection and stand **** (detail: stabilising bow)




Footrest with footrest adapter



Seat wedge, below (+ 10°), max. 15° seat tilt with additional use of rest position



Table



Belt extension (+ 20 cm / 7.9") incl. longer crotch strap



Crotch pad (large)

Product & accessories

	item code
RECARO Sport HERO Reha	880
Footrest, short* (depth: 12 cm / 4.7") (tiltable and adjustable in angle)	806
Footrest pad (10 cm / 3.9")	E11050
Footrest adapter*	821
Turning plate**** with footrest adapter**	820
Recommended options for stability:	
▪ Seatfix connection*** (for turning plate) option 1	818
▪ Stabilising bow, rear*** (for turning plate) option 2	819
▪ Stand, front (for turning plate)	817
Seat wedge, below** (+ 10°)	810
Table	885
Crotch pad (large)	E91035
Belt extension (+ 20 cm / 7.9") incl. longer crotch strap	824
Car seat protection mat	813
Car seat protection mat incl. back	814

* The adapter must be ordered to use the footrest

** Turning plate and seat wedge (below) cannot be used in combination

*** Seatfix connection and stabilising bow (rear) cannot be used in combination

**** Turning plate cannot be used when fastening the RECARO Sport HERO Reha car seat with the 3-point-belt (only ECE group I: 9 kg / 19.8" up to 18 kg / 45.7") of the vehicle. Recommended options for stability – You should choose 819 or 818 when ordering a turning plate 820!