

RECARO Monza Nova 2 Reha Order form



**Rehab Safety Car Seat
for children**

**Restraint system
for special needs „S“**

ECE R 44/04 - Group II - III

**from 15 kg / 33.1 lbs
tested up to 50 kg / 110.2 lbs**

**Order fax:
+49 (0)4761 886-19**

Status: 01.08.2021
valid until further notice

Address

Customer no.:

Company:

Street:

Postal Code:

City:

Contact person:

Phone:

E-mail:

Date:

BASIC MODEL

item code

RECARO Monza Nova 2 Seatfix Reha

9800008013000

Seat and back with integrated 11-fold height adjustable headrest; including integrated Seatfix locking arms which can be fixed to the location points of vehicles which are equipped with Isofix; stepless neck support of the headrest by an integrated air cushion; side impact protection by energy absorbing elements; removable head, lateral trunk and side supports; fixed armrests; 5-point positioning belt; crotch pad; removable abduction block; integrated speakers in the headrest for mp3-player use, integrated mp3-player pocket; RECARO Air Circulation System for optimum air circulation; ventilated foam; removable covers, can be machine-washed at 30°C (86°F)

Measurements: seat depth: 30 cm / 11.8" (43 cm / 16.9" with seat depth extension) / seat width: front 29 cm / 11.4", back 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside) / back height: 60 - 77 cm / 23.6 - 30.3"

ACCESSORIES

item code

Adapter is required to use the + accessories.

Adapter (footrest / seat depth extension / stand)

9800008212000

+ **Seatfix connection** (for swivel base)** *option 1* (if vehicle is equipped with Isofix)

9800008181000

+ **Stabilising bow, rear** (for swivel base) *option 2* (if the vehicle is NOT equipped with Isofix)

9800008190000

↓
can be combined with:

+ **Footrest, long*** (depth: 18 cm / 7.1"), tiltable and adjustable in angle

9800008051010

+ **Footrest, short*** (depth: 12 cm / 4.7"), tiltable and adjustable in angle

9800008060010

+ **Seat depth extension*** (+13 cm / 5.1")

9800008082000

+ **Stand, front** (for swivel base)** (recommended when using „Seatfix“, not necessary when using the footrest)

9800008170010

Swivel base with adapter is required to use the + accessories.

Swivel base with adapter*** (incl. adapter plate)

9800008202000

+ **Back stabilisers** (required when using swivel base)

9800008161000

+ **Seatfix connection** (for swivel base)** *option 1* (if vehicle is equipped with Isofix)

9800008181000

+ **Stabilising bow, rear** (for swivel base) *option 2* (if the vehicle is NOT equipped with Isofix)

9800008190000

↓
can be combined with:

+ **Footrest, long*** (depth: 18 cm / 7.1"), tiltable and adjustable in angle

9800008051010

+ **Footrest, short*** (depth: 12 cm / 4.7"), tiltable and adjustable in angle

9800008060010

+ **Stand, front** (for swivel base)* (recommended when using „Seatfix“, not necessary when using the footrest)

9800008170010

↪ please find further accessories on page 2

RECARO Monza Nova 2 Reha Order form

Company _____

Contact person: _____

Status: 01.08.2021
 valid until further notice

Order fax:
+49 (0)4761 886-19

ACCESSORIES	item code
<input checked="" type="checkbox"/> Abduction block	<i>included in basic model</i>
<input type="checkbox"/> Table	9800008110012
<input type="checkbox"/> Soft harness strap cover (pair)	9800008151000
<input type="checkbox"/> Chest clip	9800004230000
<input type="checkbox"/> Seat wedge, inside	9800008090011
<input type="checkbox"/> Seat wedge, below* (+10°)	9800008103011
<input type="checkbox"/> Car seat protection mat incl. back	9800008141000

Hint for more safety: With the anti-escape vest children remain belted securely in the car seat.

Anti-escape vest **9800008260000**

NOTE:

* Only footrest OR seat depth extension OR stand can be used.

** In some vehicles, the combination of seat wedge (below) and Seatfix connection cannot be used.

*** When using the swivel base, the Seatfix locking arms on the basic model may have to be removed depending on the vehicle type.