

- Head padding
- 5-point positioning harness
- Lateral trunk supports
- Crotch pad
- Lateral side supports
- Abduction block



Rehab-equipment
INCLUSIVE 

ACCESSORIES



The practical padded swivel base, the Seatfix connection for ISOFIX and the footrest give additional stability.



The seat wedge inside and below provide a relaxed resting position during the journey.



The table is a support surface and an additional impact shield that also prevents the child from attempting to unbuckle the 5-point positioning harness.



EXTRA SUPPORTIVE
 for children with special needs

Measurements & weights	RECARO Monza Nova 2 Reha
Seat depth	30 cm / 11.8" (43 cm / 16.9" with seat depth extension)
Seat width	front: 29 cm / 11.4" // rear: 24 cm / 9.4" (20 cm / 7.9" with seat wedge, inside)
Back height (seat)	60 - 77 cm / 23.6 - 30.3"
Shoulder height (child)	33 - 48 cm / 13 - 18.9"
Lower leg length (footrest mounted on swivel plate)	23 - 42 cm / 9.1 - 16.5"
Lower leg length (footrest mounted directly on the seat)	18.5 - 37.5 cm / 7.3 - 14.8"
Overall dimensions (w x h x d)	54 x 66 x 42 cm / 21.3 x 26 x 16.5"
Weight of seat approx.	7 kg / 15.4 lbs
Carrying capacity	15 - 50 kg / 33.1 - 110.2 lbs
Children height	94 - 150 cm / 37 - 59"

Products & accessories	item code
RECARO Monza Nova 2 Seatfix Reha (for Isofix)	801
Abduction block	basic configuration
Swivel base with footrest adapter** (incl. adapter plate)	820/2
Recommended options for stability:	
▪ Seatfix connection (for swivel base) option 1	818/1
▪ Stabilising bow, rear (for swivel base) option 2	819
▪ Stand, front (for swivel base)	817
Back stabilisers (pair) (required when using swivel base)	816
Footrest, long* (depth: 18 cm / 7.1") (tiltable and adjustable in angle)	805
Footrest, short* (depth: 12 cm / 4.7") (tiltable and adjustable in angle)	806
Seat depth extension* (+ 13 cm / + 5.1")	808
Seat wedge, inside	809
Seat wedge, below** (+ 10°)	810
Table	811
Anti-escape vest, size 2	823
Car seat protection mat	813
Car seat protection mat incl. back	814
Soft harness strap cover, grey (pair)	815

* Footrest and seat depth extension cannot be used in combination, unless the swivel base is used in addition.

** Swivel base and seat wedge (below) cannot be used in combination.



Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG
Walkmühlenstraße 1 · D - 27432 Bremervörde
Phone: +49 (0)4761 886-63 or -68 · Fax: +49 (0)4761 886-19
www.thomashilfen.com





RECARO Monza Nova 2 Reha
Car seat for children with special needs

“ The safety of my child is most important to me. Generally, I am a careful driver, but you never know. In case of an accident I have to be 100% sure that the car seat protects my child against injuries.

Susanne Bessel, mother


“ **RECARO** Monza Nova 2 Reha is reliably tested many times

- ✓ UN/ECE R 44/04
(United Nations Economic Commission for Europe – regulation for child restraint systems)
- ✓ DEKRA (tests for „Special Needs Restraint“)
incl. testing of all rehab accessories
- ✓ FMVSS 213 (US-American safety standard)
- ✓ Manufacturer’s guarantee



Pay attention to this seal

Thus you can be sure that the **rehab seat** successfully passed the necessary overload-test up to 50 kg / 110.2 lbs according to ECE R 44/04.



RECARO Monza Nova 2 Reha

Lightweight car seat exceeds highest safety standard



Quick and easy seat transfer
from one vehicle to another.

Extra-firm armrests
for even greater
protection in the
event of side impact.



Lightest rehab seat in its
class that you can carry
with one finger ...
Yet it carries a child up to
15 to 50 kg / 33 to 110 lbs.





Extra-stable anti-static side pads reduce crash-generated forces to the head, neck and shoulders by up to 30%.

- 30%

Stay seated!

With the anti-escape vest children remain belted securely in the car seat. Reduces danger of wriggling out of the belt.

(accessory)



Thomashilfen  Finely designed to make life easier