

NEW



EASyS Advantage S

All-round use –
concentrated on the basics


Thomashilfen

We got back to basics:

This is **NEW** ...

Well-thought-out functionality and ease of use – with the most important and significant ways of making settings, the EASyS Advantage S rehab stroller is the right choice for standard care.

This makes it the perfect companion in everyday all-round use.

Choose between 5 upholstery styles:

red/ grey	blue/ grey	green/ grey	magenta/ grey	grey/ grey

Cushions are breathable and washable

Chassis colour: black



140° resting position for children who need it



Easy seat assembly and rapid change of viewing direction – without tools:

- forward facing
- rear facing



A single release lever on the back

- seat tilt-in-space
- separates the seat from the base frame



Transportation lock to prevent unwanted unfolding, e.g. when lifting the stroller into the car

No tools needed for seat width and depth adjustments

Extremely thick cushions for great comfort.



Soft suspension absorbs vibrations and ensures a good riding pleasure without compromising steering



Convenient brake, clearly visible in red



New sportier footrest



Lever for easy adjustment of the knee angle at the pivot point

Much more storage space available under the seat:

- washable basket (inclusive)
- vent tray (optional)

EASyS Advantage S

All-round use – concentrated on the basics

Chassis & Accessories (selection)



Q Chassis:
for kindergarten, school
and at home with
comfortable height and seat
tilt-in-space adjustment



Handy accessory plate
for medical devices, e.g.
oxygen tank or suction
equipment



Grip rail and flexible
butterfly/chest harness for
a comfortable and secure
position of the upper body



Flexible lateral trunk
supports and a 4-point
pelvis belt for positioning
the upper body and the
pelvis



Handy tip assist to facilitate
curbs and other obstacles



Abduction block for
guiding the thigh and
padded h-belt for holding
the upper body



Comfortable sleeping bag,
padded inside or with
woven fur for cold days

**You will find the total
overview of all rehab
accessories in the order
form.**

Technical data	EASyS Advantage S 1	EASyS Advantage S 2
Seat depth	19 - 30 cm 7.5 - 11.8"	27 - 38 cm 10.6 - 15"
Seat width	19 - 31 cm 7.5 - 12.2"	25 - 37 cm 9.8 - 14.6"
Seat height	47 cm / 18.5"	46 cm / 18.1"
Back height	41.5 - 68 cm / 16.3 - 26.8"	55 - 74 cm / 21.7 - 29.1"
Lower leg length	17.5 - 29 cm / 6.9 - 11.4"	21 - 37 cm / 8.3 - 14.6"
Hip angle	90° to 140°	90° to 140°
Knee angle	90° to 180°	90° to 180°
Foot angle	78° to 100°	78° to 100°
Turning range of seat	180°	180°
Size of hip guide (w x d)	23 x 12 cm 9.1 x 4.7"	30 x 15 cm 11.9 x 5.9"
Size of foot rest (w x d)	28 x 18 cm 11 x 7.1"	32 x 21 cm 12.6 x 8.3"
Weight seat unit	8.5 kg / 18.7 lb	9.8 kg / 21.6 lb
max. user weight seat unit	35 kg / 77.2 lb	40 kg / 88.2 lb

Technical data	A1 chassis	A2 chassis
Seat tilt-in-space in running order	-10° to +35°	-10° to +35°
Seat tilt-in-space face-to-face	-10° to +35°	-10° to +35°
Push bar height	71.5 - 117 cm / 28.1 - 46"	71.5 - 117 cm / 28.1 - 46"
Wheel size swivel wheels (front)	18 cm / 7.1"	18 cm / 7.1"
Wheel size (rear)	25 cm / 9.8"	25 cm / 9.8"
Weight chassis	9.8 kg / 21.6 lb	10.9 kg / 24 lb
max. carrying capacity basket	5 kg / 11 lb	5 kg / 11 lb
max. carrying capacity (incl. seat)	50 kg / 110.2 lb	55 kg / 121.2 lb

Technical data stroller	EASyS Advantage S 1	EASyS Advantage S 2
Overall dimensions (l x w x h)	99.5 x 59.5 x 114 cm 39.2 x 23.4 x 44.9"	107 x 69 x 120 cm 42.1 x 27.2 x 47.2"
Folded dimensions (l x w x h)	82 x 59.5 x 64 cm 32.3 x 23.4 x 25.2"	86 x 69 x 75 cm 33.9 x 27.2 x 29.5"