

# EASyS Modular S Order form



Rehab stroller for children with short back and various head rests

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

**Order fax:**  
**+49 (0)4761 886-19**

Commission:

**Status: 28.11.2018**  
valid until further notice

COMPLETE STROLLER	item code size 1	item code size 2
-------------------	------------------	------------------

EASyS Modular S seat unit + A chassis Advantage, swivel front wheels with transit tie-down system  6208/1  6408/1






EASyS Modular S seat unit	item code size 1	item code size 2
---------------------------	------------------	------------------



EASyS Modular S seat unit (without chassis)  6210/1  6410/1

Please choose **one** pair of side paddings per seat unit / complete stroller **w/o extra charge.**

Configuration COMPLETE STROLLER / SEAT UNIT

	<b>Red side padding, pair</b>	<input type="checkbox"/> 6242	<input type="checkbox"/> 6442
	<b>Blue side padding, pair</b>	<input type="checkbox"/> 6243	<input type="checkbox"/> 6443
	<b>Green side padding, pair</b>	<input type="checkbox"/> 6244	<input type="checkbox"/> 6444
	<b>Magenta side padding, pair</b>	<input type="checkbox"/> 6241	<input type="checkbox"/> 6441
	<b>Grey side padding, pair</b>	<input type="checkbox"/> 6240	<input type="checkbox"/> 6440

Please choose **one** belt, **partly with surcharge.**


<b>2-point-pelvic belt, with front pull</b> (w/o surcharge)	<input type="checkbox"/> 8420/1	<input type="checkbox"/> 8470/1
<b>4-point-pelvic belt, with front pull</b> (with surcharge)	<input type="checkbox"/> 8421/1	<input type="checkbox"/> 8471/1

Please choose **one** headrest **with extra charge.**

<b>Headrest, shell</b> (with surcharge)	<input type="checkbox"/> 6261/1	<input type="checkbox"/> 6261/1
<b>Headrest, neck support</b> (with surcharge)	<input type="checkbox"/> 6262/1	<input type="checkbox"/> 6262/1
<b>Headrest, adjustable</b> (not approved for use as a seat in motor vehicles) (with surcharge)	<input type="checkbox"/> 6267	<input type="checkbox"/> 6267
<b>Headrest fitting, square bar with ball joint</b> (for custom-made solutions) (with surcharge)	<input type="checkbox"/> 6268	<input type="checkbox"/> 6268

**Standard features:**

**EASyS Modular S seat unit:** physiological pivots for the hip and knee joints; headrest options (selectable); breathable seat and back upholstery as well as side padding (selectable); additional foam pads for more sitting comfort; fire retardant fabrics; lateral body guides; tiltable footrest, height and angle adjustable; brackets for sun and rain cover; continuously adjustable back height, hip angle, seat depth and seat width; easily accessible adjustments without any tools needed

 accessories & chassis see page 2 & 3

Customer no.: 

Status: 28.11.2018

Order fax:

valid until further notice

+49 (0)4761 886-19

Company: Contact: 

ACCESSORIES EASyS Modular S seat unit	item code size 1	item code size 2
<b>FEET</b>		
Foot straps, black	<input type="checkbox"/> 6648/2	<input type="checkbox"/> 6648/2
Calf strap, black	<input type="checkbox"/> 6635/5	standard
<b>SEAT</b>		
Abduction block, grey	<input type="checkbox"/> 6620/6	<input type="checkbox"/> 6720/6
Grip rail, grey	<input type="checkbox"/> 6621/5	<input type="checkbox"/> 6721/5
Grip rail, grey, strong curved	<input type="checkbox"/> 6621/51	-
Tray, transparent	<input type="checkbox"/> 6628/1	<input type="checkbox"/> 6728/1
Lateral trunk support with chest belt, black	<input type="checkbox"/> 6633/2	<input type="checkbox"/> 6633/2
Lateral trunk support, grey	<input type="checkbox"/> 6634/6	<input type="checkbox"/> 6734/6
Frame padding, black	<input type="checkbox"/> 6622/6	<input type="checkbox"/> 6622/6
Side padding, red, pair	<input type="checkbox"/> 6242	<input type="checkbox"/> 6442
Side padding, blue, pair	<input type="checkbox"/> 6243	<input type="checkbox"/> 6443
Side padding, green, pair	<input type="checkbox"/> 6244	<input type="checkbox"/> 6444
Side padding, magenta, pair	<input type="checkbox"/> 6241	<input type="checkbox"/> 6441
Side padding, grey, pair	<input type="checkbox"/> 6240	<input type="checkbox"/> 6440
Footrest-tube upholstery, grey	<input type="checkbox"/> 6230	<input type="checkbox"/> 6230
<b>HARNESSES</b>		
2-point pelvic belt with front pull	<input type="checkbox"/> 8420/1	<input type="checkbox"/> 8470/1
4-point pelvic belt with front pull	<input type="checkbox"/> 8421/1	<input type="checkbox"/> 8471/1
Crotch strap (fits for 8420, 8421, 8470, 8471, 6629/5)	<input type="checkbox"/> 6341/1	<input type="checkbox"/> 6341/1
Soft pelvic harness, black	<input type="checkbox"/> 6614/6	<input type="checkbox"/> 6714/6
Butterfly / chest harness, black (S)	<input type="checkbox"/> 6616/6	-
Butterfly / chest harness, black (M)	<input type="checkbox"/> 6617/6	<input type="checkbox"/> 6617/6
Butterfly / chest harness, black (L)	-	<input type="checkbox"/> 6717/6
Butterfly / chest harness, flexible (S)	<input type="checkbox"/> 6385/3	-
Butterfly / chest harness, flexible (M)	<input type="checkbox"/> 6585/3	<input type="checkbox"/> 6585/3
Seat waistcoat, black (S)	<input type="checkbox"/> 6631/6	-
Seat waistcoat, black (M)	<input type="checkbox"/> 6632/6	<input type="checkbox"/> 6632/6
Seat waistcoat, black (L)	-	<input type="checkbox"/> 6731/6
H-belt, padded	<input type="checkbox"/> 6629/5	<input type="checkbox"/> 6629/5
<b>HEAD</b>		
Headrest, shell	<input type="checkbox"/> 6261/1	<input type="checkbox"/> 6261/1
Headrest, neck support	<input type="checkbox"/> 6262/1	<input type="checkbox"/> 6262/1
Headrest, adjustable (not approved for use as a seat in motor vehicles)	<input type="checkbox"/> 6267/1	<input type="checkbox"/> 6267/1
Headrest fitting, square bar with ball joint (for custom-made solutions)	<input type="checkbox"/> 6268	<input type="checkbox"/> 6268
<b>WEATHER PROTECTION</b>		
Sleeping bag, padded (summer), black	<input type="checkbox"/> 6619/5	<input type="checkbox"/> 6719/5
Sleeping bag, woven fur (winter), black	<input type="checkbox"/> 6618/5	<input type="checkbox"/> 6718/5
Rain shield for sun and rain cover, transparent	<input type="checkbox"/> 6739	<input type="checkbox"/> 6739
Transparent cover	<input type="checkbox"/> 6638/1	<input type="checkbox"/> 6638/1
Sun and rain cover with leg blanket	<input type="checkbox"/> 6637/6	<input type="checkbox"/> 6637/6
Sun and rain cover without leg blanket	<input type="checkbox"/> 6637/61	<input type="checkbox"/> 6637/61
Umbrella, grey	<input type="checkbox"/> 6594	<input type="checkbox"/> 6594

# EASyS Modular S Order form (chassis & accessories)

**Status: 28.11.2018**
**Order fax:**

valid until further notice

**+49 (0)4761 886-19**

Company: \_\_\_\_\_

Contact: \_\_\_\_\_

**A chassis Advantage for outdoor use**
**item code size 1**
**item code size 2**

**A chassis Advantage**, swivel front wheels with transit tie down system

 6202/3

 6402/3

**Standard features:**
**A chassis Advantage:** lightweight and stable aluminum chassis, transit tie-down system; foldable; one-hand adjustment of seat tilt-in-space; face-to-face option; rear frame suspension; height adjustable push bar with biocompatible soft grip; convenient foot brake; no-flat EVA-wheels; swivel front wheels incl. swivel lock; removable and washable basket; folding lock

**ACCESSORIES for A chassis Advantage**
**item code size 1**
**item code size 2**
**Vent tray**
 6294/1

 6494/1

**Pneumatic wheels** (upgrade kit)

 6647/1

 6647/1

**Tip assist, anthracite**
 6750/2

 6750/2

**Q chassis for indoor use**

**Q chassis** (indoor use)

 6666/6

 6666/6

**Adapter for Q chassis** (necessary in combination with EASyS Modular S seat unit)

 6295/1

 6295/1

**Standard features:**
**Q chassis:** comfortable indoor base, easy height and seat tilt-in-space adjustment with gas spring (foot engaged), 4 swivelling castors (2 lockable)

**ACCESSORIES for Q chassis**
**item code size 1**
**item code size 2**
**Push bar, black**
 6660/1

 6660/1