

EASyS Modular Order form



Rehab stroller for children
with safety seat

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

Order fax:
+49 (0)4761 886-19

Status: 26.01.2018
valid until further notice

COMPLETE STROLLER	item code size 1	item code size 2
-------------------	------------------	------------------

EASyS Modular seat unit + A chassis, swivel front wheels with swivel lock

6209

6409

EASyS Modular seat unit + A chassis, swivel front wheels with transit tie down system

6208

6408

EASyS Modular seat unit	item code size 1	item code size 2
-------------------------	------------------	------------------



EASyS Modular seat unit (without chassis)

6210

6410

Please choose **one** pair of side paddings and **one** headrest per seat unit / complete stroller **w/o** extra charge.

SIDE PADDING

Configuration COMPLETE STROLLER / SEAT UNIT



Berry side padding, pair

6263

6463



Classic side padding, pair

6264

6464

HEADREST

Headrest, standard black

6260

6460

Headrest, shell black

6261

6261

Headrest, neck support black

6262

6262

BELT

2-point-pelvic belt, with front pull (w/o surcharge)

8420

8470

4-point-pelvic belt, with front pull (with surcharge)

8421

8471

Standard features:

EASyS Modular seat unit: BSS® (Balance Safety System) - safety seat; physiological pivots for the hip and knee joint; adapter for different headrests; headrest (freely selectable); adjustable arm supports (height, width and depth); breathable seat and back upholstery with reflective safety-seams; side padding (freely selectable); fire retardant fabrics; lateral body guides; tiltable footrest, height and angle adjustable; EASyS click-system for an easy seat installation / removal; brackets for sun and rain cover; continuously adjustable back height, hip angle, seat depth and seat width; easily accessible adjustments using an Allen key

A chassis	item code size 1	item code size 2
-----------	------------------	------------------

A chassis, swivel front wheels with swivel lock

6201

6401

A chassis, fixed front wheels

6200

6400

A chassis, swivel front wheels with transit tie down system

6202

6402

Standard features:

A chassis: lightweight and stable aluminum chassis, foldable; tilt-in-space lever at push bar; lockable seat suspension; continuously adjustable seat tilt-in-space; height adjustable push bar with biocompatible soft grip; foot brake; no-flat EVA-wheels; fixed front wheels or swivel front wheels incl. swivel lock; removable basket

Indoor bases	item code size 1 + 2
--------------	----------------------

Q hi/low base standard

6666

Z standard underframe for indoor use

6691

 accessories see page 2

EASyS Modular Order form (accessories)

Status: 26.01.2018
Order fax:

valid until further notice

+49 (0)4761 886-19

Company: _____

Contact: _____

ACCESSORIES EASyS Modular seat unit	item code size 1	item code size 2
FEET		
Foot straps, black	<input type="checkbox"/> 6648/1	<input type="checkbox"/> 6648/1
Calf strap, black	<input type="checkbox"/> 6635/5	standard
SEAT		
Abduction block, grey	<input type="checkbox"/> 6620/6	<input type="checkbox"/> 6720/6
Grip rail, grey	<input type="checkbox"/> 6621/5	<input type="checkbox"/> 6721/5
Grip rail, grey, strong curved	<input type="checkbox"/> 6621/51	-
Tray, transparent	<input type="checkbox"/> 6628/1	<input type="checkbox"/> 6728/1
Lateral trunk support with chest belt, black	<input type="checkbox"/> 6633/2	<input type="checkbox"/> 6633/2
Lateral trunk support, grey	<input type="checkbox"/> 6634/6	<input type="checkbox"/> 6734/6
Frame padding, black	<input type="checkbox"/> 6622/6	<input type="checkbox"/> 6622/6
Side padding, berry, pair	<input type="checkbox"/> 6263	<input type="checkbox"/> 6463
Side padding, classic, pair	<input type="checkbox"/> 6264	<input type="checkbox"/> 6464
HARNESSES		
2-point pelvic belt with front pull	<input type="checkbox"/> 8420/1	<input type="checkbox"/> 8470/1
4-point pelvic belt with front pull	<input type="checkbox"/> 8421/1	<input type="checkbox"/> 8471/1
Soft pelvic harness, black	<input type="checkbox"/> 6614/6	<input type="checkbox"/> 6714/6
Butterfly / chest harness, black (S)	<input type="checkbox"/> 6616/6	-
Butterfly / chest harness, black (M)	<input type="checkbox"/> 6617/6	<input type="checkbox"/> 6617/6
Butterfly / chest harness, black (L)	-	<input type="checkbox"/> 6717/6
Butterfly / chest harness, flexible (S)	<input type="checkbox"/> 6385/3	-
Butterfly / chest harness, flexible (M)	<input type="checkbox"/> 6585/3	<input type="checkbox"/> 6585/3
Seat waistcoat, black (S)	<input type="checkbox"/> 6631/6	-
Seat waistcoat, black (M)	<input type="checkbox"/> 6632/6	<input type="checkbox"/> 6632/6
Seat waistcoat, black (L)	-	<input type="checkbox"/> 6731/6
H-belt, padded	<input type="checkbox"/> 6629/5	<input type="checkbox"/> 6629/5
HEAD		
Headrest, standard, black	<input type="checkbox"/> 6260	<input type="checkbox"/> 6460
Headrest, shell, black	<input type="checkbox"/> 6261	<input type="checkbox"/> 6261
Headrest, neck support, black	<input type="checkbox"/> 6262	<input type="checkbox"/> 6262
WEATHER PROTECTION		
Sleeping bag, padded (summer)	<input type="checkbox"/> 6619/5	<input type="checkbox"/> 6719/5
Sleeping bag, woven fur (winter)	<input type="checkbox"/> 6618/5	<input type="checkbox"/> 6718/5
Rain shield for sun and rain cover, transparent	<input type="checkbox"/> 6739	<input type="checkbox"/> 6739
Transparent cover	<input type="checkbox"/> 6638/1	<input type="checkbox"/> 6638/1
Sun and rain cover with leg blanket	<input type="checkbox"/> 6637/6	<input type="checkbox"/> 6637/6
Sun and rain cover without leg blanket	<input type="checkbox"/> 6637/61	<input type="checkbox"/> 6637/61
Umbrella	<input type="checkbox"/> 6594	<input type="checkbox"/> 6594

ACCESSORIES for A chassis	item code size 1 + 2
Tip assist, anthracite	<input type="checkbox"/> 6750/1
Transit-tie-down system	<input type="checkbox"/> E6783/9

ACCESSORIES for Q chassis	item code size 1 + 2
Push bar for Q hi/low base	<input type="checkbox"/> 6660