



Rehab stroller for children

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

**Order fax:**  
+49 (0)4761 886-19

**Status: 26.05.2021**  
valid until further notice

BASIC MODEL	item code size 1	item code size 2
-------------	------------------	------------------



Swiftly stroller

9800068065000

**Standard features Swiftly:** lightweight and foldable chassis with integrated seat unit; crash eyelets, H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover

**Measurements:**

- seat depth: 22 cm / 8.7" (17 cm / 6.9" with seat minimizer) - 28,5 cm / 11.2"
- seat width: 34 cm / 13.4" (23 cm / 9.1" with seat minimizer)
- back height: 62 cm / 24.4"
- min. shoulder height: 37 cm / 14.6" (26 cm / 10.2" with seat minimizer)



Swiftly 2 stroller

9800068080010

**Standard features Swiftly 2:** lightweight and foldable chassis with integrated seat unit; crash eyelets, H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover; folding lock

**Measurements:**

- seat depth: 34 cm / 13.4" (24.5 cm / 9.6" with seat minimizer) - 40 cm / 15.7"
- seat width: 39.5 cm / 15.6" (29.5 cm / 11.6" with seat minimizer)
- back height: 70 cm / 27.6"
- min. shoulder height: 47 cm / 18.5" (30 cm / 11.8" with seat minimizer)

 accessories see page 2

Customer no.:

Status: 26.05.2021

valid until further notice

Company:

Contact:

ACCESSORIES	item code size 1	item code size 2
<b>FEET / LEGS</b>		
Foot straps black (pair)	<input type="checkbox"/> 9800066482000	<input type="checkbox"/> 9800066482000
Footpad (Seat minimizer)	<input type="checkbox"/> 5005100239005	-
Basket	<input type="checkbox"/> 9800068130000	<input type="checkbox"/> 9800068350000
Tip assist, anthracite	<input type="checkbox"/> 9800068291000	<input type="checkbox"/> 9800068291000
<b>SEAT</b>		
Abduction block, grey (not possible in combination with crotch strap)	<input type="checkbox"/> 9800066206000	<input type="checkbox"/> 9800067206000
Grip rail, grey	<input type="checkbox"/> 9800068331000	<input type="checkbox"/> 9800067215000
Tray, transparent	<input type="checkbox"/> 9800068220000	<input type="checkbox"/> 9800067281010
Lateral trunk support, grey	<input type="checkbox"/> 9800066346000	<input type="checkbox"/> 9800067346000
Lateral trunk support with chest belt	<input type="checkbox"/> 9800066332000	<input type="checkbox"/> 9800066332000
Lateral seat upholstery, red	<input type="checkbox"/> 5005100228000	<input type="checkbox"/> 9700442011000
Frame padding black	<input type="checkbox"/> 9800066226000	<input type="checkbox"/> 9800066226000
Frame padding, red-black	<input type="checkbox"/> 9800066227000	<input type="checkbox"/> 9800066227000
Seat minimizer, set black	<input type="checkbox"/> 9800068340000	<input type="checkbox"/> 9800068360000
Seat minimizer, set red	<input type="checkbox"/> 9800068341000	-
Seat width reducer, black (seat width: 29.5 cm / 11.6")	-	<input type="checkbox"/> 9800068410000
Seat width reducer, red (seat width: 29.5 cm / 11.6")	-	<input type="checkbox"/> 9800068411000
Backrest upholstery, black (Seat minimizer)	<input type="checkbox"/> 5005100219005	-
Backrest upholstery, red (Seat minimizer)	<input type="checkbox"/> 5005100226000	-
Backrest upholstery, black (Seat minimizer)	<input type="checkbox"/> 5005100229005	-
Backrest upholstery, red (Seat minimizer)	<input type="checkbox"/> 5005100227000	-
<b>HARNESSES</b>		
2-point pelvic harness with front-pull	<input type="checkbox"/> 9800084201000	<input type="checkbox"/> 9800084701000
4-point pelvic harness with front-pull	-	<input type="checkbox"/> 9800084711000
Crotch strap (can be used only in combination with „pelvic harness“ or „H-belt“)	<input type="checkbox"/> 9800063411000	<input type="checkbox"/> 9800063411000
Soft pelvic harness, black	<input type="checkbox"/> 9800066146000	<input type="checkbox"/> 9800067146000
Butterfly/chest harness, black (S)	<input type="checkbox"/> 9800066166000	-
Butterfly/chest harness, black (M)	<input type="checkbox"/> 9800066176000	<input type="checkbox"/> 9800066176000
Butterfly/chest harness, black (L)	-	<input type="checkbox"/> 9800067176000
Butterfly/chest harness, flexible (S)	<input type="checkbox"/> 9800063853000	-
Butterfly/chest harness, flexible (M)	<input type="checkbox"/> 9800065853000	<input type="checkbox"/> 9800065853000
Seat waistcoat, black (S)	<input type="checkbox"/> 9800066315100	-
Seat waistcoat, black (M)	<input type="checkbox"/> 9800066325100	<input type="checkbox"/> 9800066325100
Seat waistcoat, black (L)	-	<input type="checkbox"/> 9800067315100
<b>HEAD</b>		
Head pillow, black	<input type="checkbox"/> 9800068240000	<input type="checkbox"/> 9800068370000
Head pillow, red	-	<input type="checkbox"/> 9800068371000
Head pillow, black (Seat minimizer)	<input type="checkbox"/> 5005100209005	-
Head pillow, red (Seat minimizer)	<input type="checkbox"/> 5005100225000	<input type="checkbox"/> 9800068420000
<b>WEATHER PROTECTION</b>		
Sleeping bag, padded (body contour), summer	<input type="checkbox"/> 9800066195000	<input type="checkbox"/> 9800067195000
Sleeping bag, woven fur (body contour), winter	<input type="checkbox"/> 9800066185000	<input type="checkbox"/> 9800067185000
Rain cover, transparent, with leg blanket	<input type="checkbox"/> 9800066381000	<input type="checkbox"/> 9800066381000
Sun and rain cover with leg blanket	<input type="checkbox"/> 9800066376000	<input type="checkbox"/> 9800066376000
Sun and rain cover without leg blanket	<input type="checkbox"/> 9800066376100	<input type="checkbox"/> 9800066376100
Rain shield for sun and rain cover, transparent (can be used only in combination with „sun and rain cover“)	<input type="checkbox"/> 9800067390000	<input type="checkbox"/> 9800067390000
Umbrella	<input type="checkbox"/> 9800065940000	<input type="checkbox"/> 9800065940000
<b>CHASSIS</b>		
Transit tie down system	<input type="checkbox"/> 9800068300000	standard
Pneumatic wheels (upgrade kit)	<input type="checkbox"/> 9800066471000	<input type="checkbox"/> 9800066710000