



Rehab stroller for children

Address

Customer no.:

Company: _____

Street: _____

ZIP Code / Place: _____

Phone: _____

Date: _____

Contact: _____

Order fax:
+49 (0)4761 886-19

Status: 26.01.2018
valid until further notice

BASIC MODEL	item code size 1	item code size 2
-------------	------------------	------------------



- Swiftly 1 stroller standard version**
- Swiftly 1 stroller with swivel lock front wheels**
- Swiftly 1 stroller with transit tie-down-kit**

- 6806/3**
- 6807/3**
- 6806/31**

Standard features Swiftly 1: lightweight and foldable chassis with integrated seat unit; H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; continuously adjustable seat depth; swivel front wheels; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover



Swiftly 2 stroller with swivel lock front wheels & transit tie-down-kit

6808

Standard features Swiftly 2: lightweight and foldable chassis with integrated seat unit; H-belt; fire retardant fabrics; breathable upholstery; hip-angle adjustable backrest; swivel front wheels; no-flat EVA-wheels; height adjustable and tiltable footrest; foot brake; height adjustable push bar; bracket for sun and rain cover; folding lock

 accessories see page 2

Customer no.:

Status: 26.01.2018

valid until further notice

Company:

Contact:

ACCESSORIES	item code size 1	item code size 2
FEET / LEGS		
Foot straps black (pair)	<input type="checkbox"/> 6648/1	<input type="checkbox"/> 6648/2
Basket	<input type="checkbox"/> 6813	<input type="checkbox"/> 6835
Tip assist, anthracite	<input type="checkbox"/> 6829/1	<input type="checkbox"/> 6829/1
SEAT		
Abduction block, grey	<input type="checkbox"/> 6620/6	<input type="checkbox"/> 6720/6
Grip rail, grey	<input type="checkbox"/> 6833/1	<input type="checkbox"/> 6721/5
Tray, transparent	<input type="checkbox"/> 6822	<input type="checkbox"/> 6728/1
Seat minimizer, set	<input type="checkbox"/> 6834	<input type="checkbox"/> 6836
Lateral trunk support	<input type="checkbox"/> 6634/6	<input type="checkbox"/> 6734/6
Lateral trunk support with chest belt	<input type="checkbox"/> 6633/2	<input type="checkbox"/> 6633/2
Frame padding black	<input type="checkbox"/> 6622/6	<input type="checkbox"/> 6622/6
HARNESSES		
2-point pelvic harness with front-pull	<input type="checkbox"/> 8420/1	<input type="checkbox"/> 8470/1
4-point pelvic harness with front-pull	–	<input type="checkbox"/> 8471/1
Crotch strap (fits to 8420, 8470, H-belt)	<input type="checkbox"/> 6341	<input type="checkbox"/> 6341/1
Soft pelvic harness, black	<input type="checkbox"/> 6614/6	<input type="checkbox"/> 6714/6
Butterfly/chest harness, black (S)	<input type="checkbox"/> 6616/6	–
Butterfly/chest harness, black (M)	<input type="checkbox"/> 6617/6	<input type="checkbox"/> 6617/6
Butterfly/chest harness, black (L)	–	<input type="checkbox"/> 6717/6
Butterfly/chest harness, flexible (S)	<input type="checkbox"/> 6385/2	–
Butterfly/chest harness, flexible (M)	<input type="checkbox"/> 6585/2	<input type="checkbox"/> 6585/3
Seat waistcoat, black (S)	<input type="checkbox"/> 6631/6	–
Seat waistcoat, black (M)	<input type="checkbox"/> 6632/6	<input type="checkbox"/> 6632/6
Seat waistcoat, black (L)	–	<input type="checkbox"/> 6731/6
HEAD		
Head pillow	<input type="checkbox"/> 6824	<input type="checkbox"/> 6837
WEATHER PROTECTION		
Sleeping bag, padded (body contour), summer	<input type="checkbox"/> 6619/5	<input type="checkbox"/> 6719/5
Sleeping bag, woven fur (body contour), winter	<input type="checkbox"/> 6618/5	<input type="checkbox"/> 6718/5
Rain cape, grey	<input type="checkbox"/> 6340/4	<input type="checkbox"/> 6540/4
Rain cover, transparent	<input type="checkbox"/> 6638/1	<input type="checkbox"/> 6638/1
Sun and rain cover with leg blanket	<input type="checkbox"/> 6637/6	<input type="checkbox"/> 6637/6
Sun and rain cover without leg blanket	<input type="checkbox"/> 6637/61	<input type="checkbox"/> 6637/61
Rain shield for sun and rain cover, transparent (can be used only in combination with 6637/6 or 6637/61)	<input type="checkbox"/> 6739	<input type="checkbox"/> 6739
Umbrella	<input type="checkbox"/> 6594	<input type="checkbox"/> 6594
CHASSIS		
Transit tie down system	<input type="checkbox"/> 6830	standard
Pneumatic wheels (upgrade kit)	<input type="checkbox"/> 6647	<input type="checkbox"/> 6671