

ThevoTwist Order form



**Rehab stroller
with mobility seat
for children**

Address

Customer no.:

Company:

Street:

ZIP Code / Place:

Phone:

Date:

Contact:

**Order fax:
+49 (0)4761 886-19**

Status: 05.02.2019
valid until further notice

COMPLETE STROLLER	item code size 1	item code size 2
-------------------	------------------	------------------

ThevoTwist seat unit + A chassis, swivel front wheels with swivel lock	<input type="checkbox"/> 6213 (red) <input type="checkbox"/> 6212 (blue)	<input type="checkbox"/> 6413 (red) <input type="checkbox"/> 6412 (blue)
ThevoTwist seat unit + A chassis, swivel front wheels with transit tie-down system	<input type="checkbox"/> 6207 (red) <input type="checkbox"/> 6206 (blue)	<input type="checkbox"/> 6407 (red) <input type="checkbox"/> 6406 (blue)



ThevoTwist seat unit (ThevoSeat)	item code size 1	item code size 2
----------------------------------	------------------	------------------



ThevoTwist seat unit (ThevoSeat)
(without chassis)

6688 6788

Please choose upholstery color **without surcharge** and headrest and belt, **partly with surcharge**.

Configuration COMPLETE STROLLER / SEAT UNIT	COLOR	 <input type="checkbox"/> red <input type="checkbox"/> red  <input type="checkbox"/> blue <input type="checkbox"/> blue
	HEADREST	<input type="checkbox"/> 6387/1 <input type="checkbox"/> 6587/1 <input type="checkbox"/> 6387 <input type="checkbox"/> 6587 <input type="checkbox"/> 6599/1 <input type="checkbox"/> 6599/1 <input type="checkbox"/> 6599 <input type="checkbox"/> 6599 <input type="checkbox"/> 6598/1 <input type="checkbox"/> 6598/1 <input type="checkbox"/> 6598 <input type="checkbox"/> 6598
	BELT	<input type="checkbox"/> 8420 <input type="checkbox"/> 8470 <input type="checkbox"/> 8421 <input type="checkbox"/> 8471
	Standard features:	ThevoTwist seat unit: BSS (Balance Safety System) - safety seat; flexible wing elements and wing suspensions in the back; physiological pivots for the hip and knee joint; hip angle adjustment with gas spring; breathable seat and back upholstery with reflective safety-seams; fire retardant fabrics; brackets for sun and rain cover; belt bow; femoral support; tiltable footrest, height and angle adjustable; EASyS click-system for an easy seat installation/removal; easily accessible adjustments using an Allen key

A chassis	item code size 1	item code size 2
-----------	------------------	------------------

A chassis, swivel front wheels with swivel lock	<input type="checkbox"/> 6201/2	<input type="checkbox"/> 6401/2
A chassis, fixed front wheels	<input type="checkbox"/> 6200/2	<input type="checkbox"/> 6400/2
A chassis, swivel front wheels with transit tie down system	<input type="checkbox"/> 6202	<input type="checkbox"/> 6402

Standard features:

A chassis: lightweight and stable aluminum chassis, foldable; tilt-in-space lever at push bar; lockable seat suspension; continuously adjustable seat tilt-in-space; height adjustable push bar with biocompatible soft grip; foot brake; no-flat EVA-wheels; fixed front wheels or swivel front wheels incl. swivel lock; removable basket

Indoor bases	item code size 1 + 2
--------------	----------------------

Q hi/low base standard	<input type="checkbox"/> 6666
Z standard underframe for indoor use	<input type="checkbox"/> 6691

 accessories see page 2

ThevoTwist Order form (accessories)

Customer no.:

Status: 05.02.2019**Order fax:**

valid until further notice

+49 (0)4761 886-19

Company:

Contact:

ACCESSORIES ThevoTwist seat unit	item code size 1	item code size 2
FOOTREST		
Foot straps, black	<input type="checkbox"/> 6648/2	<input type="checkbox"/> 6648/2
Calf strap, grey	<input type="checkbox"/> 6381	-
SEAT		
Abduction block, grey	<input type="checkbox"/> 6620/6	<input type="checkbox"/> 6720/6
Side plate, high	<input type="checkbox"/> 6386	<input type="checkbox"/> 6586
Frame padding, black	<input type="checkbox"/> 6622/6	<input type="checkbox"/> 6622/6
Grip rail, grey	<input type="checkbox"/> 6372	<input type="checkbox"/> 6572
Grip rail, grey, strong curved	<input type="checkbox"/> 6372/1	-
Tray, transparent	<input type="checkbox"/> 6628/1	<input type="checkbox"/> 6728/1
Lateral trunk support plane	<input type="checkbox"/> 6377/1	<input type="checkbox"/> 6577/1
Lateral trunk support with chest belt	<input type="checkbox"/> 6596	<input type="checkbox"/> 6596
Impulse sensors (lateral trunk support MiS)	<input type="checkbox"/> 6383	<input type="checkbox"/> 6583
Moldable aluminum lath	<input type="checkbox"/> 6384	<input type="checkbox"/> 6584
Upholstery for impulse sensors (has to be ordered in combination with moldable aluminium lath)	<input type="checkbox"/> 500512042	<input type="checkbox"/> 500512043
HARNESSES		
2-point pelvic belt with front pull	<input type="checkbox"/> 8420/1	<input type="checkbox"/> 8470/1
4-point pelvic belt with front pull	<input type="checkbox"/> 8421/1	<input type="checkbox"/> 8471/1
Soft pelvic harness, black	<input type="checkbox"/> 6614/6	<input type="checkbox"/> 6714/6
Seat waistcoat, black (S)	<input type="checkbox"/> 6631/51	<input type="checkbox"/> -
Seat waistcoat, black (M)	<input type="checkbox"/> 6632/51	<input type="checkbox"/> 6632/51
Seat waistcoat, black (L)	-	<input type="checkbox"/> 6731/51
Butterfly/chest harness, black (S)	<input type="checkbox"/> 6616/6	-
Butterfly/chest harness, black (M)	<input type="checkbox"/> 6617/6	<input type="checkbox"/> 6617/6
Butterfly/chest harness, black (L)	-	<input type="checkbox"/> 6717/6
Butterfly/chest harness, flexible (S)	<input type="checkbox"/> 6385/3	<input type="checkbox"/> -
Butterfly/chest harness, flexible (M)	<input type="checkbox"/> 6585/3	<input type="checkbox"/> 6585/3
H-belt, padded	<input type="checkbox"/> 6629/5	<input type="checkbox"/> 6629/5
Belt bow (has to be ordered in combination with butterfly/chest harness + h-belt)	<input type="checkbox"/> 6371/2	<input type="checkbox"/> 6571/2
HEADREST		
Head rest, standard, red	<input type="checkbox"/> 6387/1	<input type="checkbox"/> 6587/1
Head rest, standard, blue	<input type="checkbox"/> 6387	<input type="checkbox"/> 6587
Head rest, neck support, red	<input type="checkbox"/> 6598/1	<input type="checkbox"/> 6598/1
Head rest, neck support, blue	<input type="checkbox"/> 6598	<input type="checkbox"/> 6598
Head rest, shell, red	<input type="checkbox"/> 6599/1	<input type="checkbox"/> 6599/1
Head rest, shell, blue	<input type="checkbox"/> 6599	<input type="checkbox"/> 6599
Headrest height reducer	<input type="checkbox"/> 6569	<input type="checkbox"/> 6569
WEATHER PROTECTION		
Sleeping bag, padded (summer)	<input type="checkbox"/> 6373	<input type="checkbox"/> 6573
Sleeping bag, woven fur (winter)	<input type="checkbox"/> 6374	<input type="checkbox"/> 6574
Rain cape, grey	<input type="checkbox"/> 6340/4	<input type="checkbox"/> 6540/4
Rain shield for sun and rain cover, transparent	<input type="checkbox"/> 6739	<input type="checkbox"/> 6739
Transparent cover	<input type="checkbox"/> 6638/1	<input type="checkbox"/> 6638/1
Sun and rain cover with leg blanket, black	<input type="checkbox"/> 6637/6	<input type="checkbox"/> 6637/6
Sun and rain cover without leg blanket, black	<input type="checkbox"/> 6637/61	<input type="checkbox"/> 6637/61
Umbrella	<input type="checkbox"/> 6594	<input type="checkbox"/> 6594

Note:
The belt bow has to be ordered in combination with the butterfly/chest harness + h-belt, otherwise the mounting of the belts is not possible.

ACCESSORIES for A chassis	item code size 1 + 2
Pneumatic wheels (upgrade kit)	<input type="checkbox"/> 6647
Tip assist, anthracite	<input type="checkbox"/> 6750/1
Transit-tie-down system	<input type="checkbox"/> E67839

ACCESSORIES for Q chassis	item code size 1 + 2
Push bar for Q hi/low base	<input type="checkbox"/> 6660