

Recommended for
walking therapy!

trio

■ walker ■ ■



Now with
adjustable width!


Thomashilfen

3 advantages of the walker

- + puncture safe „Flexel“-tyres
- + individual height adjustment
- + arm rests adjustable in width

Further advantages:

- high stability and tilt safety, tested according to DIN EN ISO 11199-3-2005
- adjustable without tools
- ideal for assembly of purpose built items (e.g. pelottes)
- height adjustable armpit supports (optional)
- lockable hand brakes
- small turn radius
- walking practise inside and outside
- available in 2 different versions



Height adjustable and tilting armpit supports (optional).



Comfortable arm rest, height and width adjustable.



Adjustable handles and lockable hand brake included.

Product description	Item code	length	width	height (arm rest)	wheel-Ø	max. weight capacity
Walker trio B (+ width adjustment)	2045	90 cm 35.4"	71 - 77 cm 28 - 30.3"	90 - 130 cm 35.4 - 51.2"	200 x 50 mm 7.9 x 2"	150 kg 331 lb.
Walker trio BA (+ width adjustment + armpit support)	2040	90 cm 35.4"	71 - 77 cm 28 - 30.3"	90 - 130 cm 35.4 - 51.2"	200 x 50 mm 7.9 x 2"	150 kg 331 lb.